

# Isotopic Resolution of Chromatographically Separated IdeS Subunits Using the X500B QTOF System

Fan Zhang<sup>2</sup>, Sean McCarthy<sup>1</sup>  
<sup>1</sup>SCIEX, MA, USA, <sup>2</sup>SCIEX, CA USA

Analysis of protein subunits using high resolution mass spectrometry (HRMS) is commonly performed during the development of biotherapeutics. In comparison to intact mass analysis, subunit level approaches offer greater ability to chromatographically separate subunits and improved sensitivity to detect and quantify a wide range of post translational modifications.

Recently for subunit analysis there has been a drive towards collection of mass spectrometric data which has been isotopically resolved to provide greater confidence in peak identification. The X500B QTOF is specifically designed for generation of high quality data across a number of biopharmaceutical workflows including intact mass, subunit and peptide analysis for a wide range of therapeutic modalities.

Presented in this technical note is the use of the SCIEX X500B QTOF System for confirmation of antibody subunits generated by digestion with IdeS. This molecule is used as an example of the resolution possible on this platform for proteins or their subunits within a similar mass range. Mass spectrometric and data processing parameters to achieve high quality data are presented and discussed.



## Key Feature of X500B QTOF Solution

- High resolution mass spectrometer for a wide range of biopharmaceutical applications
- Isotopically resolved subunit analysis by HRMS
- Compact benchtop footprint reduces laboratory space requirements
- Easy to use hardware and software accessible for a wide range of users




Figure 1. Comparison of Theoretical and Experimental Reconstructed Light Chain Spectra from IdeS Digested NIST Antibody.

## Methods

### Sample Preparation:

IdeS reduced NISTmAb was first digested with IdeS protease (Cat V7511, Promega) following the protocol from Promega. Briefly, 1 unit of IdeS protease was incubated with 1 µg of IgG at 37°C for 1h in PBS. The resultant sample was further reduced with TCEP using the above protocol. 5 µl samples were subjected to LC-MS analysis.

### Chromatography:

Separation was accomplished using a Shimadzu ExionLC™ System fitted with an Agilent PLRP-S column (2.1mm X 50mm, 300Å, 5µm) at 80°C using the gradient shown in Table 1. Mobile phase A was 0.1% formic acid in water and mobile phase B was 0.1% formic acid in acetonitrile.

**Table 1. LC Gradient Conditions**

Time (min)	%A	%B	Flow Rate ml/min
Initial	75	25	0.25
3.0	75	25	0.25
9.0	10	90	0.25
11.4	10	90	0.25
11.5	75	25	0.25
17.0	75	25	0.25

### Mass Spectrometry:

A SCIEX X500B Mass Spectrometer with a Turbo V™ Source was used for data acquisition. Data was acquired using TOF-MS mode with intact protein mode (IPM) turned off. MS instrument conditions are listed in Table 2.

**Table 2. MS Parameters.**

Parameter	Setting
Scan Mode	Positive
GS1	50
GS2	50
Curtain Gas	35
Temperature	400°C
Ion Spray Voltage	5000 V
Time Bins to Sum	6
Accumulation Time (ms)	0.5 sec
TOF Start Mas (Da)	400
TOF Stop Mas (Da)	3000
Declustering Potential	150.0
Collision Energy	10

### Data Processing:

Data were processed using SCIEX OS software 1.4.0.18067 and Bio Tool Kit software 1.0. Reconstruction parameters are shown in Table 3

**Table 3. Reconstruction Parameters.**

Parameter	Setting
Start Mass	Fc 25000 Da LC 22600 Da Fd 25200 Da
Stop Mass	Fc 26000 Da LC 23600 Da Fd 26200 Da
Step Mass	0.05 Da
Input Spectrum Isotope Resolution	30000

## Results and Discussion

IdeS digestion and reduction of a therapeutic monoclonal antibody results in three roughly equivalent mass fragments which correspond to the Fc/2, light chain, and Fd fragments of the monoclonal antibody. This digestion strategy has gained significant attention since it was first reported for use on therapeutic monoclonal antibodies.<sup>1</sup> In particular, IdeS digestion serves to reduce overall sample complexity compared to intact mass measurement and reduces the complexity of sample preparation and data interpretation frequently associated with peptide mapping experiments.

Decreasing sample using and IdeS digest enables greater interrogation of the protein subunit when compared to intact mass analysis. As each fragment generated is around 25 kDa it is possible to generate isotopically resolved spectral data for each component. As shown in Figures 2-4, under the conditions described, isotopically resolved spectra for each subunit were obtained with major relevant features observed compared to previous reports<sup>2</sup>




Figure 2. Reconstructed Spectrum for Fc/2 Fragment of NISTmAb Standard Digested with IdeS.





Figure 3. Reconstructed Spectrum for Light Chain of NISTmAb Standard Digested with IdeS and Reduction.



**Figure 4. Reconstructed Spectrum for Fd of NISTmAb Standard Digested with IdeS and Reduction.**

To verify the quality of our data, we generated a theoretical spectrum based on the elemental composition of the light chain of the NISTmAb standard. The theoretical spectrum was modeled to generate a continuous spectrum with the same resolution used for reconstruction of the experimental data, ca. 30,000. As shown in Figure 1, excellent agreement of experimental (blue trace) and theoretical (green trace) was achieved.

Finally, consistent with previous reports<sup>2</sup>, the monoisotopic mass for any species observed is not clearly evident as its abundance is too low in intensity to accurately define. In this case, we verified the mass accuracy of our measurement by comparison of theoretical isotopic masses to those observed experimentally for the light chain. Shown in Table 4 are the theoretical and observed reconstructed masses for each isotope clearly visible in the data. The average error across the light chain was 3.46 PPM highlighting the high mass accuracy with which subunit data can be characterized with the X500B QTOF system.

## Conclusions

- Isotopically resolved IdeS subunit analysis on the X500B QTOF system provides accurate assessment of post translational modification with high mass accuracy.
- Strong agreement of theoretical and experimental spectral data for subunits is achieved.
- Presented data is highly consistent with previous reports.

**Table 4. Theoretical and Observed masses for light chain with corresponding PPM values.**

Theoretical	Observed	PPM
23118.318	23118.222	-4.13
23119.320	23119.222	-4.25
23120.323	23120.226	-4.19
23121.325	23121.237	-3.82
23122.328	23122.238	-3.89
23123.330	23123.25	-3.48
23124.333	23124.256	-3.33
23125.335	23125.257	-3.39
23126.338	23126.255	-3.58
23127.340	23127.257	-3.60
23128.343	23128.264	-3.40
23129.345	23129.267	-3.37
23130.347	23130.273	-3.22
23131.350	23131.277	-3.14
23132.352	23132.281	-3.07
23133.354	23133.292	-2.70
23134.357	23134.291	-2.84
23135.359	23135.294	-2.81
	<b>AVERAGE</b>	<b>3.46</b>

## References

1. An, Y., Zhang, Y., Mueller, H.-M., Shameem, M., Chen, X.,  
A New Tool for Monoclonal Antibody  
Analysis. *mAbs*. 2014;6(4) :879-893. doi:  
10.4161/mabs.28762.
2. John E. Schiel, Darryl L. Davis, Oleg V. Borisov (Eds.). Vol  
1201, 2015 State-of-the-Art and Emerging Technologies for  
Therapeutic Monoclonal Antibody Characterization  
Volume2. Biopharmaceutical Characterization: The  
NISTmAb Case Study. Vol 1201, American Chemical  
Society.

AB Sciex is doing business as SCIEX.

© 2018 AB Sciex. For Research Use Only. Not for use in diagnostic procedures. The trademarks mentioned herein are the property of AB Sciex Pte. Ltd. or their respective owners. AB SCIEX™ is being used under license.

Document number: RUO-MKT-02-8241-A



### Headquarters

500 Old Connecticut Path | Framingham, MA 01701 USA  
Phone 508-383-7700  
[sciex.com](http://sciex.com)

### International Sales

For our office locations please call the division  
headquarters or refer to our website at  
[sciex.com/offices](http://sciex.com/offices)