

Forensic Identification and Quantification Workflows Delivered on a Revolutionary Designed QTOF and SCIEX OS Software

Igniting your routine forensic testing with the new SCIEX X500R QTOF System

Xiang He¹ and Adrian Taylor,²

¹SCIEX, 1201 Radio Rd, Redwood City, CA 94065, USA; ²SCIEX, 71 Four Valley Drive, Concord, Ontario, L4K 4V8 Canada.

Overview

Quadrupole Time-of-Flight (QTOF) mass spectrometry is becoming the desired technology for sensitive and selective screening workflows in a forensic toxicological setting. The technology overcomes many challenges faced when using traditional techniques and more significantly captures all information about the sample in one injection to allow for retrospectively mining the data. Using the accurate mass and mass resolution information from both TOF-MS and TOF-MS/MS acquired data allows for simultaneous highly specific targeted quantitation and non-targeted screening. Here we describe a new benchtop QTOF system with revolutionary N geometry TOF designed flight path and new, intuitive software for easy adoption of accurate mass technology to forensic testing. We demonstrate that the new hardware and software combined allow a high level of confidence for compound identification and quantification from urine samples in one seamless workflow.

Introduction

Liquid Chromatography coupled to Tandem Mass Spectrometry (LC-MS/MS) is a widely used analytical tool for the screening of compounds and metabolites. Triple quadrupole based mass analyzers operated in Multiple Reaction Monitoring (MRM) mode have become the preferred method to routinely deliver highly selective and sensitive quantitative results, but are limited to targeted screening only.

With an increasing demand for retrospective and non-targeted analyses of forensic toxicological samples, high resolution, accurate mass and full scan mass analyzers are gaining popularity. The adoption of the technology has been restricted by more complicated to use and more expensive instrumentation compared to their nominal mass counterparts. Here we introduce a revolutionary new Quadrupole Time-of-Flight (QTOF) mass spectrometer that contains advances in engineering design to bring the high performance TOF-MS and TOF-MS/MS capabilities into a compact benchtop platform.



Figure 1: The SCIEX ExionLC™ AC HPLC system (left), the SCIEX X500R QTOF System (middle) and SCIEX OS Software (right).

The SCIEX X500R QTOF mass spectrometer is part of a complete workflow from the fully integrated SCIEX ExionLC™ Systems to the freshly designed SCIEX OS software; a new user interface for simultaneous identification and quantification workflows (Figure 1.)

SCIEX X500R QTOF System

The new benchtop SCIEX X500R QTOF System with revolutionary N geometry TOF designed flight path has been engineered for simplicity, service accessibility and minimized footprint. N TOF geometry, versus V geometry, gives the same effective flight path length for ions and therefore resolution, but in a smaller overall foot print. This has been accomplished with an extra mirror in the TOF chamber without a loss in transmission (Figure 2). To maintain stable mass accuracy the system uses a simple heated TOF vacuum chamber design. This consists of 6 discreet heater drones maintaining a constant TOF chamber temperature, insulating against ambient temperature changes (Figure 2).




Figure 2. SCIEX X500R QTOF System and Technology Advances

The system has been designed to maximize robustness and uptime

- Integrated Calibrant Delivery System and Turbo V™ Source with TwinSpray probe (Figure 3), allows seamless mass accuracy auto-calibrations during long runs.
- Service Accessibility
 - Easy QJet® and Turbo pump access for fast and efficient maintenance, increasing system uptime
 - Segmented TOF vacuum chamber – allows easy access to detector while protecting sensitive accelerator.



Figure 3. Integrated Calibrant Delivery System and Turbo V™ Source with TwinSpray probe

Figure 4 shows the mass accuracy stability of the SCIEX X500R QTOF System when analyzing multiple urine samples, spiked with various concentrations of analytes, without auto-calibration, over a ten hour period. The majority of compounds are shown to be within a 1 ppm mass accuracy over this time period.




Figure 4. Mass Accuracy Stability of the SCIEX X500R QTOF System in the Analysis of Urine Samples

Figure 5 shows the resolution for both TOF-MS and TOF-MS/MS masses sampled over a seven day time period on a SCIEX X500R QTOF System.




Figure 5. Resolution of the SCIEX X500R QTOF System Over a Week's Period for Selected m/z ; both TOF-MS and TOF-MS/MS

Figure 6 shows a representative linear dynamic range of the SCIEX X500R QTOF System showing 4 orders for the Asenapine compound.

SCIEX OS Software

SCIEX OS Software is a single software platform for LC and MS control, data processing as well as reporting.



Figure 6. Linearity of the SCIEX X500R QTOF System shown for Asenapine (0.5 ng/mL to 1000 ng/mL)

The SCIEX OS software is intuitive and logical, segregated into Acquisition, Processing and Management work spaces (Figure 7). In the Acquisition work space there are separate method editors for the LC and MS parameters as well as batch creation and queue panes. The Processing allows for simultaneous identification and quantification. The Management workspace allows the adjustment of hardware, software and user settings.




Figure 7. Home Page of SCIEX OS Software. Single Software Platform for LC/MS Control, Data Processing and Reporting.

Acquisition

The SCIEX OS software has a simplified step by step acquisition method setup with only relevant parameters being visible. Figure 8 shows the setup for an Information Dependent Acquisition method for the analysis of small molecules and the intuitive steps that are taken to input the MS parameter values.

For a quick instrument status check, the Manual Tune guides the user through the steps to perform a quick review of the system performance, perform an auto-calibration and report out the test result prior to running a batch (Figure 9).




Figure 8. SCIEX OS Software MS Acquisition Method




Figure 9. SCIEX OS Software MS Tune Allows for Quick Instrument Status Check via Simple Step by Step Instructions

Building a batch is assisted by the smart grid design allowing copy/paste, fill down, auto increment and import/export. Figure 10 shows the batch editor and the link to the auto-calibration setup.

Once the batch has been submitted to the queue the Auto-Cal samples are inserted as shown in the Queue Manager in Figure 11. The SCIEX OS software allows for detailed instrument status including monitoring and recording of LC pressure traces as well as direct control of the individual components of the system (Figure 11).




Figure 10. SCIEX OS Software Batch Editor and Setup for Auto-Calibration




Figure 12. SCIEX OS Software Allows the Simultaneous Review of Qualitative and Quantitative Results

The SCIEX OS Software allows the user to filter the results to only show compounds that pass acceptance criteria and are detected with user defined confidence (Figure 13)




Figure 11. SCIEX OS Software Queue Manager with Inserted Auto-Cal Samples and Detailed Instrument Status Panel

Processing- Analytics

Once a results table is generated, quantitative and qualitative results can be reviewed in the same panel (Figure 12). A Traffic light system indicates the confidence of the identification based on accurate mass, retention time, isotopic pattern and library matching. Compounds calculated to be above the cutoff concentration in unknown samples are flagged. In the same workspace the peak integration, spectra and calibration lines can be displayed.



Figure 13. SCIEX OS Software Filtering Criteria

Finally results can be reported out using the *Create Report* functionality (Figure 14)




Figure 14. SCIEX OS Software Report Generation

Acquisition Workflows on the SCIEX X500R QTOF System with SCIEX OS Software

Information Dependent Acquisition

Information Dependent Acquisition (IDA) is a non-targeted data acquisition (Figure 15). It allows for TOF-MS quantification and provides high confidence in screening with MS/MS information that uses high selectivity through unit Q1 resolution. IDA-MS/MS provides the most interference-free fragmentation information.



Figure 15. Information Dependent Acquisition

When creating an IDA acquisition, MS and MS/MS settings are all contained in a single User Interface. Figure 16 shows the parameters used in the IDA experiments described in this technical note. In this example, one TOF-MS survey scan and up to 16 dependent TOF-MS/MS scans are triggered from the survey scan, in each data cycle.

Figure 16. SCIEX OS Software Information Dependent Acquisition Method Editor

Due to the high scanning speed (up to 100 Hz for single collision energy) on SCIEX X500R QTOF systems, almost all potential

compound targets in the samples can be confirmed with confident MS/MS library matching.

IDA-MS/MS is a non-targeted data acquisition method and the user needs to define the maximum number of candidates in each data cycle. More intense ions take higher priority within any data cycle, so for less abundant species especially in complex sample matrices, the associated MS/MS information might be missed. Therefore, an unbiased MS/MS data acquisition approach that collects MS/MS information for everything at all times (MS/MS^{All}) will solve this potential concern.

SWATH[®] Acquisition

SWATH[®] acquisition (Figure 17) is non-targeted and provides MS/MS information for everything in the sample, all the time. Each scan cycle in SWATH[®] Acquisition starts with a TOF-MS experiment. The acquisition approach therefore allows for screening and quantification from both TOF-MS and TOF-MS/MS acquired data.

Most of the existing MS/MS^{All} techniques collect MS and MS/MS information for all ions in an alternating fashion, i.e. MS scan of all precursor ions, followed by MS/MS scan of the fragments of all precursor ions. Without precursor ion selection, such approaches suffer from insufficient sensitivity, selectivity and narrower linear range compared to IDA-MS/MS.




Figure 17. SWATH[®] Acquisition

SWATH[®] acquisition uses either a fixed or a variable Q1 isolation window, as part of a TOF-MS/MS experiment, which is stepped across the mass range of interest. Figure 18 shows the SWATH[®] acquisition method editor in the SCIEX OS Software, with the example of 16 looped TOF-MS/MS experiments, each with a different (variable) Q1 isolation window, that are required to cover the mass range of interest (120 to 500 *m/z*).




Figure 18. SCIEX OS Software SWATH® Acquisition Method Editor

By varying the Q1 isolation window for each TOF-MS/MS experiment we are able to separate compounds with similar mass into different SWATH® Acquisition windows so that we minimize the amount of convolution (multiple precursor ions generating common fragment ions at the same time) in each TOF-MS/MS experiment (Figure 19).




Figure 19. Constructing Variable SWATH® Acquisition Window Sizes for each Looped TOF-MS/MS Experiment to Minimize Convolution in the SCIEX OS Software

MRM^{HR}

MRM^{HR} (High Resolution Multiple Reaction Monitoring) is a targeted data acquisition for quantification purposes and can be unscheduled or scheduled. Compound dependent parameters can be optimized for each MRM^{HR}.

MRM (SCIEX Triple Quad® or QTRAP®)



MRM^{HR} (QTOF: TripleTOF® or X500R QTOF)




Figure 20. Comparison of MRM^{HR} with traditional (unit resolution) MRM

To help transition familiarity of MRM performed on a triple quadrupole to MRM^{HR} performed on the SCIEX X500R QTOF system, the SCIEX OS Software has a unique way of building the MRM^{HR} method to have the look and feel of performing traditional MRM experiments by allowing the input of the precursor ion mass (MRM Q1 equivalent mass) and accurate fragment mass (MRM Q3 equivalent nominal mass) (Figure 21). These transitions can easily be imported from the SCIEX high resolution 1700 compound MS/MS forensic spectral library to include up to 5 transitions per compound.




Figure 21. SCIEX OS Software Scheduled MRM^{HR} Method Editor, Fragment Ion Mass ± 10 m/z

The quantification method is then generated automatically from the acquisition method (Figure 22).

Row	IS	Group	Name	Precursor Ma...	Fragment Mas (Da)	XIC Width (Da)	Retention Time (min)	IS--	Experiment Index
1		6-MAM	6-MAM	328.154	211.0747	0.02	2.04		2 TOF MSMS of 328.2 (40 - 500)
2		7-Aminoclonazepam	7-Aminoclonazepam	286.074	222.1025	0.02	2.67		3 TOF MSMS of 286.1 (40 - 500)
3		7-Hydroxymirtazapine	7-Hydroxymirtazapine	435.223	190.0684	0.02	2.96		4 TOF MSMS of 435.2 (40 - 500)
4		Acetyl Fentanyl	Acetyl Fentanyl	323.212	105.07	0.02	3.09		5 TOF MSMS of 323.2 (40 - 500)
5		Alpha-Hydroxyprazosin	Alpha-Hydroxyprazosin	325.085	297.0665	0.02	4.19		6 TOF MSMS of 325.1 (40 - 500)
6		Alpha-Hydroxymidazolam	Alpha-Hydroxymidazolam	342.08	203.0377	0.02	4.51		7 TOF MSMS of 342.1 (40 - 500)
7		Alpha-Hydroxytriazolam	Alpha-Hydroxytriazolam	359.046	331.0272	0.02	3.94		8 TOF MSMS of 359.0 (40 - 500)
8		Alpha-PPP	Alpha-PPP	294.138	105.0699	0.02	2.04		8 TOF MSMS of 294.1 (40 - 500)
9		Alpha-PVP	Alpha-PVP	212.17	161.0963	0.02	2.61		10 TOF MSMS of 232.2 (40 - 500)
10		Alprazolam	Alprazolam	389.09	281.073	0.02	4.49		11 TOF MSMS of 309.1 (40 - 500)
11		AM-2201 4-OH pentyl	AM-2201 4-OH pentyl	376.171	155.0492	0.02	6.07		12 TOF MSMS of 376.2 (40 - 500)
12		Amisulpride	Amisulpride	276.19	117.0702	0.02	4.53		13 TOF MSMS of 278.2 (40 - 500)
13		Amphetamine	Amphetamine	136.112	91.0553	0.02	1.63		14 TOF MSMS of 136.1 (40 - 500)
14		Benazepine	Benazepine	290.139	168.1021	0.02	2.36		15 TOF MSMS of 290.1 (40 - 500)
15		Buphedrone	Buphedrone	178.123	131.07	0.02	1.98		16 TOF MSMS of 178.1 (40 - 500)
16		Buprenorphine	Buprenorphine	448.311	414.2636	0.02	3.67		17 TOF MSMS of 448.3 (40 - 500)
17		Carisoprodol	Carisoprodol	243.181	55.0565	0.02	1.68		18 TOF MSMS of 243.2 (40 - 500)
18		Clomipramine	Clomipramine	335.162	86.0959	0.02	5.26		19 TOF MSMS of 335.2 (40 - 500)
19		Codone	Codone	360.159	215.1109	0.02	1.87		20 TOF MSMS of 360.2 (40 - 500)
20		Cotinine	Cotinine	177.102	80.0496	0.02	1.83		21 TOF MSMS of 177.1 (40 - 500)
21		Cyclobenzaprine	Cyclobenzaprine	276.175	215.0878	0.02	4.28		22 TOF MSMS of 276.2 (40 - 500)
22		Desalkylflurazepam	Desalkylflurazepam	289.054	140.0294	0.02	4.42		23 TOF MSMS of 289.1 (40 - 500)
23		Desipramine	Desipramine	247.186	72.0823	0.02	4.28		24 TOF MSMS of 247.2 (40 - 500)
24		Desmethyloxepin	Desmethyloxepin	246.154	107.0493	0.02	3.60		25 TOF MSMS of 246.2 (40 - 500)
25		Dextromethorphan	Dextromethorphan	272.201	215.1438	0.02	3.44		26 TOF MSMS of 272.2 (40 - 500)
26		Diazepam	Diazepam	285.079	154.0424	0.02	5.43		27 TOF MSMS of 285.1 (40 - 500)

Figure 22. Automatically generating the SCIEX OS Software MRM^{HR} Quantification Method from the SCIEX OS Software MRM^{HR} Acquisition Method

Alternatively, if the fragment masses are not known at the time of the acquisition method creation, then the traditional MRM^{HR} setup is still achievable by inputting the TOF start and stop masses (Figure 23).

Group name	Precursor L.	TOF sta...	TOF sta...	Accumul...	Decks...	Collision...	Retention L...	Retention time tolerance (+/- sec)
1	6-MAM	328.15	40.00000	500.00000	0.0100	60	33	2.04
2	7-Aminoclonazepam	286.07	40.00000	500.00000	0.0100	60	32	2.67
3	7-Hydroxymirtazapine	435.22	40.00000	500.00000	0.0100	60	40	2.96
4	Acetyl Fentanyl	323.21	40.00000	500.00000	0.0100	60	48	3.09
5	Alpha-Hydroxyprazosin	325.09	40.00000	500.00000	0.0100	60	31	4.18
6	Alpha-Hydroxymidazolam	342.08	40.00000	500.00000	0.0100	60	34	4.51
7	Alpha-Hydroxytriazolam	359.05	40.00000	500.00000	0.0100	60	35	3.96
8	Alpha-PPP	294.14	40.00000	500.00000	0.0100	60	35	2.04
9	Alpha-PVP	212.17	40.00000	500.00000	0.0100	60	20	2.61
10	Alprazolam	389.09	40.00000	500.00000	0.0100	60	35	4.49
11	AM-2201 4-OH pentyl	376.17	40.00000	500.00000	0.0100	60	35	6.07
12	Amisulpride	276.19	40.00000	500.00000	0.0100	60	31	4.53
13	Amphetamine	136.11	40.00000	500.00000	0.0100	60	31	1.63
14	Benazepine	290.14	40.00000	500.00000	0.0100	60	31	2.36
15	Buphedrone	178.12	40.00000	500.00000	0.0100	60	31	1.98
16	Buprenorphine	448.31	40.00000	500.00000	0.0100	60	31	3.67
17	Carisoprodol	243.18	40.00000	500.00000	0.0100	60	31	1.68
18	Clomipramine	335.16	40.00000	500.00000	0.0100	60	31	5.26
19	Codone	360.16	40.00000	500.00000	0.0100	60	31	1.87
20	Cotinine	177.10	40.00000	500.00000	0.0100	60	31	1.83
21	Cyclobenzaprine	276.17	40.00000	500.00000	0.0100	60	31	4.28
22	Desalkylflurazepam	289.05	40.00000	500.00000	0.0100	60	31	4.42
23	Desipramine	247.19	40.00000	500.00000	0.0100	60	31	4.28
24	Desmethyloxepin	246.15	40.00000	500.00000	0.0100	60	31	3.60
25	Dextromethorphan	272.20	40.00000	500.00000	0.0100	60	31	3.44
26	Diazepam	285.08	40.00000	500.00000	0.0100	60	31	5.43
84	Zolpidem	308.18	40.00000	500.00000	0.0100	60	45	3.05
85	THC-COOH	345.21	40.00000	500.00000	0.0100	60	30	6.34

Figure 23. Scheduled MRM^{HR} Method Editor, MS/MS Full Scan

Materials and Methods

Compound list and spiking solutions

Table 1 lists the 93 compounds plus internal standards. All were procured from Cerilliant Corporation (Round Rock, TX). Two spiking solutions in methanol were prepared: one for analytes (SA) and the other for internal standards (SIS). Concentrations of all the analytes in the spiking solution SA are listed in Table 1.

Compounds in black font are in the regular panel (72 analytes) and the ones in blue font are the additional 21 analytes in the extended panel (93 analytes). Internal standards are shown in grey background.

Calibrator preparation

Blank human urine was spiked with solution SA to prepare calibrators. Four levels of calibrators were prepared. Actual concentrations varied for each compound, however the concentration ratio between these calibrators was always (in descending order): 20:6:2:1. For instance, the four different concentrations (in descending order) for fentanyl in the calibrators were: 20, 6, 2 and 1 ng/mL, while those of gabapentin were: 1000, 300, 100 and 50 ng/mL.

Sample preparation

- 100 μ L urine sample was mixed with 25 μ L IMCS Rapid Hydrolysis Buffer, 20 μ L IMCSzyme and 10 μ L SIS. Both IMCS Rapid Hydrolysis Buffer and IMCSzyme were acquired from IMCS (Columbia, SC). Hydrolysis time was typically between 30 and 60 min at 55°C.
- After hydrolysis was complete, 0.2 mL methanol and 0.625 mL water were added to the mixture.
- The mixture was then centrifuged at 21,000 g for 10 min.
- The supernatant was transferred to a glass vial with insert for analysis by LC-MS/MS.

Liquid Chromatography

Liquid Chromatography analysis was performed on the SCIEX ExionLCTM AC HPLC system at 30°C. Separation was achieved using a Phenomenex Kinetex Phenyl-Hexyl column (50 \times 2.1 mm, 2.6 μ m, 00B-4495-E0), with a Phenomenex SecurityGuard ULTRA UHPLC Phenyl (AJ0-8774) and ULTRA holder (AJ0-9000). Mobile phase A (MPA) was ammonium formate in water. Mobile phase B (MPB) was formic acid in methanol. The LC flow rate was 0.5 mL/min and the LC run-times investigated were 8.0 and 2.0 minutes. Injection volume was 10 μ L.

Table 1: List of analytes and internal standards, and their concentrations in spiking solution (for preparation of calibrators)

Compounds	(ng/mL)	Compounds	(ng/mL)	Compounds	(ng/mL)	Compounds	(ng/mL)
6-MAM	1000	Gabapentin	10000	Naloxone	5000	Pentobarbital	10000
7-Aminoclonazepam	5000	Hydrocodone	5000	Naltrexone	5000	Secobarbital	10000
7-Hydroxymitragynine	1000	Hydromorphone	5000	N-desmethyltapentadol	5000	THC-COOH	2000
Acetyl Fentanyl	200	Imipramine	5000	Norbuprenorphine	2000	6-MAM-d3	
Alpha-Hydroxyalprazolam	5000	JWH 122 5-OH pentyl	1000	Norcodeine	5000	Amphetamine-d5	
Alpha-Hydroxymidazolam	5000	JWH 19 6-OH hexyl	1000	Nordiazepam	5000	Benzoyllecgonine-d3	
Alpha-Hydroxytriazolam	5000	JWH 210 5-OH-pentyl	1000	Norfentanyl	200	Buprenorphine-d4	
Alpha-PPP	1000	JWH-018 4-OH pentyl	1000	Norhydrocodone	5000	Carisoprodol-d7	
Alpha-PVP	1000	JWH-018 pentanoic acid	1000	Normeperidine	5000	Codeine-d6	
Alprazolam	5000	JWH-073 3-OH butyl	1000	Noroxycodone	5000	Fentanyl-d5	
AM-2201 4-OH pentyl	1000	JWH-073-butanoic acid	1000	Norpropoxyphene	10000	Hydrocodone-d6	
Amitriptyline	5000	JWH-250-N-4-OH pentyl	1000	Nortriptyline	5000	Hydromorphone-d6	
Amphetamine	10000	JWH-073-butanoic acid	1000	O-Desmethyltramadol	5000	JWH 018 4-OH pentyl-d5	
Benzoyllecgonine	5000	JWH-250-N-4-OH pentyl	1000	Oxazepam	5000	JWH 019 6-OH hexyl-d5	
Buphedrone	1000	Lorazepam	5000	Oxycodone	5000	MDPV-d8	
Buprenorphine	2000	MDA	10000	Oxymorphone	5000	Meperidine-d4	
Carisoprodol	10000	MDEA	10000	PCP	2500	Mephedrone-d3	
Clomipramine	5000	MDMA	10000	Pregabalin	10000	Meprobamate-d7	
Codeine	5000	MDPV	1000	Propoxyphene	10000	Methadone-d3	
Cotinine	5000	Meperidine	5000	Protriptyline	5000	Methamphetamine-d5	
Cyclobenzaprine	5000	Mephedrone	1000	RCS4-4-OH-pentyl	1000	Methylone-d3	
Desalkylflurazepam	5000	Meprobamate	10000	Ritalinic Acid	5000	Mitragynine-d3	
Desipramine	5000	Methadone	10000	Sufentanil	200	Morphine-d6	
Desmethyldoxepin	5000	Methamphetamine	10000	Tapentadol	5000	Nordiazepam-d5	
Dextromethorphan	5000	Methedrone	1000	Temazepam	5000	Nortriptyline-d3	
Diazepam	5000	Methylone	1000	Tramadol	5000	Oxycodone-d6	
Dihydrocodeine	5000	Methylphenidate	5000	Zolpidem	5000	Oxymorphone-d3	
Doxepin	5000	Midazolam	5000	Amobarbital/pentobarbital	10000	THC-COOH-d3	
EDDP	10000	Mitragynine	1000	Butabarbital	10000	Butalbital-d5	
Fentanyl	200	Morphine	5000	Butalbital	10000	Secobarbital-d5	

Grey background: IS

SCIEX OS Software Processing

Identification and Quantification Results

Defining the retention time and accurate precursor and fragment mass for each analyte is performed first (Figure 24) followed by setting up the library searching parameters.

Row	IT	Group	Name	Chemical Formula	Adduct	Precursor Mass	Fragment Mass (Da)	XIC Width (Da)	Retention Time (min)	IS N
1	25		2,3-dimethoxy-4...	C13H21NO25	[M+H] ⁺	256.13658	107.063	0.02	1.02	
2	25		2,3-dimethoxy-4...	C13H21NO25	[M+H] ⁺	256.13658	107.063	0.02	6.90	
3	25		2,3-dimethoxy-4...	C13H21NO25	[M+H] ⁺	256.13658	124.06	0.02	6.90	
4	25		2,3-dimethoxy-4...	C13H21NO25	[M+H] ⁺	256.13658	182.04	0.02	6.90	
5	25		2,3-dimethoxy-4...	C13H21NO25	[M+H] ⁺	256.13658	167.02	0.02	6.90	
6	25		2,3-dimethoxy-4...	C13H21NO25	[M+H] ⁺	256.13658	134.01	0.02	1.74	
7	25		2,3-dimethoxy-4...	C13H21NO25	[M+H] ⁺	256.13658	166.08	0.02	6.90	
8	20-8		20-8-FLY 1	C12H14BN02	[M+H] ⁺	284.02807		0.02	0.66	
9	20-8		20-8-FLY 2	C12H14BN02	[M+H] ⁺	284.02807	188.0829	0.02	3.83	
10	20-8		20-8-FLY 3	C12H14BN02	[M+H] ⁺	284.02807	267	0.02	8.55	
11	20-8		20-8-FLY 4	C12H14BN02	[M+H] ⁺	284.02807	172.06	0.02	6.95	
12	20-8		20-8-FLY 5	C12H14BN02	[M+H] ⁺	284.02807	145.04	0.02	1.65	
13	3-4		3-4-dimethylmeth...	C12H17N02	[M+H] ⁺	192.13829		0.02	1.60	
14	3-4		3-4-dimethylmeth...	C12H17N02	[M+H] ⁺	192.13829	159.1042	0.02	4.42	
15	3-4		3-4-dimethylmeth...	C12H17N02	[M+H] ⁺	192.13829	144.08	0.02	4.83	
16	3-4		3-4-dimethylmeth...	C12H17N02	[M+H] ⁺	192.13829	145.06	0.02	4.60	
17	3-4		3-4-dimethylmeth...	C12H17N02	[M+H] ⁺	192.13829	174.11	0.02	4.43	
18	3-ites		3-desmethylprodi...	C13H21NO2	[M+H] ⁺	248.16451		0.02	4.72	
19	3-ites		3-desmethylprodi...	C13H21NO2	[M+H] ⁺	248.16451	174.1177	0.02	4.07	
20	3-ites		3-desmethylprodi...	C13H21NO2	[M+H] ⁺	248.16451	20.07	0.02	0.91	

Figure 24. Defining the Retention Time, Accurate mass of Precursor and Fragment Ions

Defining the qualifying components includes setting accuracy tolerance levels for calibrants and controls as well as flagging integration discrepancies. Qualifying definitions also includes defining the identification criteria and setting the confidence levels at which mass error, error in retention time, isotope pattern and library matching scores are deemed an acceptable difference, marginal difference or unacceptable difference (Figure 25).

Define a qualifying row:

- Integration acceptance: Any
- Accuracy: Any
- Calculated concentration: Any
- Show rows that: Qualify

Maximum tolerance for accuracy:

- Standards at Lower Limit of Quantitation (LLOQ): +/- 20.0 %
- Standards: +/- 15.0 %
- Quality Controls (QC): +/- 15.0 %

Define a qualifying row:

- Ion ratio: [check]
- Mass error: [check]
- RT: [check]
- Isotope: [check]
- Library: [check]
- C_nH_m Formula: [check]
- Show rows that: Qualify

Qualitative Rule	Acceptable Difference	Marginal Difference	Unacceptable Difference	Combined Score Weight (%)
Mass Error (ppm)	< 5	< 10	>= 10	20
Error in Retention Time	< 5	< 10	>= 10	10
% Difference Isotope Ratio	< 20	< 40	>= 40	10
Library Hit Score	> 60	> 30	<= 30	60

Figure 25. Defining the Identification and Quantification Qualifying Components in the SCIEX OS Software

Results and Discussion

As part of evaluating the new SCIEX X500R QTOF to perform simultaneous identification and quantification of compounds from forensically related samples routinely, we investigated two LC gradients. We evaluated each methods capabilities to elute all analytes throughout the entire gradient as evenly as possible in

order to maximize triggering IDA MS/MS for all components, reduce the MRM^{HR} concurrency for quality of data (*Scheduled* MRM^{HR}), resolve isobaric species and alleviate ion suppression caused by co-elution of excessive number of analytes. Figure 26 shows the Extracted Ion Chromatograms (XICs) for the 8.0 minute run and Figure 27 show the XICs for the 2.0 minute run.




Figure 26. Extracted Ion Chromatograms for Analytes from a Urine Analysis using an 8.0 Minute LC Runtime




Figure 27. Extracted Ion Chromatograms for Analytes from a Urine Analysis using an 2.0 Minute LC Runtime

Information Dependent Acquisition

With the ability to provide the most interference free fragmentation information for library searching in a non-targeted acquisition, the IDA workflow provides the highest confidence screening using MS/MS information. Figure 28 shows the multiple screening criteria that are used for identification purposes in the SCIEX OS Software's easy to understand user interface.

Component Name	Act. Conc.	Exp. RT	Height	Ret. Time	Precursor Mass	Mass Error...	RT Conf...	Isotope Conf...	Library Conf...	Found At Mass	Mass Error (ppm)	Retent. Time Er...	Library Hit	Library Score	Comb. Score	Isotope Ratio...
6-MAM	N/A	2.04	118100	2.05	328.1541		✓	✓	✓	328.1541	1.0	0.3	6-O-Monacetilmorphine	94.5	12.9	4.3
7-Aminoclonazepam	N/A	2.67	371155	2.68	286.0742		✓	✓	✓	286.0742	2.3	0.3	7-Aminoclonazepam	95.4	90.9	5.8
7-Hydroxymibogalin	N/A	2.90	23047	2.97	415.2217		✓	✓	✓	415.2217	0.8	0.4	7-Hydroxymibogalin	67.0	75.5	10.7
Acetyl Fentanyl	N/A	3.01	43934	3.01	323.2118		✓	✓	✓	323.2118	1.0	0.1	Acetyl Fentanyl	94.0	94.5	0.8
Alpha-Hydroxypropyl	N/A	4.19	100433	4.19	325.0811		✓	✓	✓	325.0811	1.0	0.0	Alpha-Hydroxypropyl	96.3	95.1	2.2
Alpha-Hydroxymid...	N/A	4.51	150403	4.51	342.0804		✓	✓	✓	342.0804	1.2	0.1	Alpha-Hydroxymidazo...	98.4	96.3	1.1
Alpha-Hydroxytraz...	N/A	3.95	85889	3.95	359.0461		✓	✓	✓	359.0461	1.0	0.0	alpha-Hydroxytriazolam	91.5	92.2	2.6
Alpha-PPP	N/A	2.04	155541	2.03	204.1363		✓	✓	✓	204.1363	1.7	0.4	Alpha-PPP	90.4	89.9	1.6
Alpha-PVP	N/A	2.61	217884	2.61	232.1696		✓	✓	✓	232.1696	1.2	0.0	Alpha-PVP	99.2	96.6	1.6
Alprazolam	N/A	4.49	393119	4.49	309.0962		✓	✓	✓	309.0962	1.8	0.1	Alprazolam	92.4	89.7	8.3
Ara-2021 4-OH Me...	N/A	6.07	121905	6.07	376.1707		✓	✓	✓	376.1709	0.4	0.1	4-Ara-2021 10-(4-Hydroxy...	100.0	99.1	5.7
Amphetamine	N/A	4.53	64880	4.52	278.1903		✓	✓	✓	278.1907	1.4	0.3	Amphetamine	98.1	95.5	1.4
Amphetamine	N/A	1.62	28842	1.64	136.1121		✓	✓	✓	136.1122	1.0	1.4	Amphetamine	100.0	91.9	10.5
Benzoylgonine	N/A	2.36	86068	2.37	290.1387		✓	✓	✓	290.1382	1.8	0.2	Benzoylgonine	99.4	95.7	0.7
Buphedrone	N/A	1.97	107591	1.98	178.1226		✓	✓	✓	178.1228	0.7	0.4	Buphedrone	99.1	97.2	1.2
Buprenorphine	N/A	3.68	111768	3.68	468.3108		✓	✓	✓	468.3112	0.7	0.1	Buprenorphine	96.3	96.0	0.9
Carisoprodol	N/A	3.66	55995	3.68	261.1809		✓	✓	✓	261.1811	1.0	0.4	Carisoprodol	96.6	95.3	1.1




Figure 28 Screening and Identification Results from an IDA Experiment

The importance of acquiring quality MS/MS data for identification purposes, and not to solely rely on the accurate mass of the precursor ion, is demonstrated in Figures 29, 30 and 31. Each figure demonstrates how, by acquiring MS/MS data, we can distinguish between structural isobaric compounds. In each example shown, isobaric compounds are barely chromatographically separated and so the presence of either or both the compounds cannot be identified by either accurate mass of the precursor ions or confidently by retention time if there is any drift in retention of the compounds. The highest confidence is gained through library MS/MS comparisons.



Figure 29. High Confidence Identification of Naloxone and 6-MAM Isobaric Compounds Gained through Library MS/MS Comparisons




Figure 30. High Confidence Identification of Buprenorphine and Mephedrone Isobaric Compounds Gained through Library MS/MS Comparisons




Figure 31. High Confidence Identification of Methylphenidate and Normeperidine Isobaric Compounds Gained through Library MS/MS Comparisons

Figures 32 and 33 show selected compound examples of XICs from the TOF-MS information acquired as part of the IDA workflow. This information can be used for quantification purposes.




Figure 32. XICs of α -PVP in Urine from TOF-MS information (Urine was diluted 10-fold; 10 μ L injection)




Figure 33. XICs of Sufentanil in Urine from TOF-MS information (Urine was diluted 10-fold; 10 μ L injection)

Figure 34 shows representative calibration curves obtained from the IDA experiment.




Figure 34. Representative Calibration Curves for Selected Compounds Showing that the TOF-MS information can be used for Quantification in an IDA Workflow

SWATH® Acquisition Results

SWATH® Acquired data can be processed in a similar way to processing IDA data for screening purposes. Again this uses multiple criteria for confidence in identification; most importantly using MS/MS library matching. Figure 35 shows a result of this from the 8.0 minute LC run which resulted in a high true positive rate of 98%.




Figure 35. Processed SWATH® Acquired Data using Multiple Identification Criteria; including MS/MS Library Matching

With traditional IDA-MS/MS, quantitation can only be performed from TOF-MS mode and not from the *in situ* sporadic TOF-MS/MS data points. In contrast, due to the continual and looped MS/MS scan function, quantification from fragment ions is achievable from SWATH® acquisition. Better selectivity from the fragment ion information (Figure 36) relative to parent ion information, allows more sensitive detection in MS/MS mode of lower concentration species in complex matrices.




Figure 36. Gains in Selectivity with the Ability to Extract Out a Specific Fragment Ion From Variable Window SWATH® Acquired Data Compared to Extracted Accurate Mass of the Precursor Ion

Figure 37 shows identification and quantification results for a synthetic drug obtained from SWATH® Acquisition using the 8.0 min LC run time. This compound was not in the original targeted list but retrospective interrogation of the data from this unknown

sample allowed for its identification without having to re-inject the sample again.



Figure 37. Identification and Quantification Results for n-Ethylcathinone Ephedrine Metabolite Compound Analysed by SWATH® Acquisition

The n-ethylcathinone ephedrine metabolite compound was identified based on unique fragment ions and their ratios as well as a library searching match (Figure 38). In a SWATH® acquisition experiment, not only can confirmation of the presence of compounds be made through MS/MS library matching and ion ratio calculations but because of the ability to extract out many unique fragment ions from the SWATH® acquired MS/MS data we can also determine the concentration based on quantification of either or both the precursor and fragment ions depending on which has less interferences.




Figure 38. Extraction of Unique Fragment Ions From SWATH® Acquisition and Using Both Ion Ratio and Library Matching to Confirm Presence of n-Ethylcathinone Ephedrine metabolite in an Unknown Urine Sample

When investigating using a 2.0 minute LC run time as part of the SWATH® acquisition, we were able to accomplish good quantification results. Figure 39 shows representative calibration curves obtained from the ultra-fast screening experiment.




Figure 39. Representative Calibration Curves Generated from the SWATH® Acquisition using a 2.0 minute LC Runtime (n=3)

Sensitivity examples are shown in Figures 40, 41 and 42 for selected compounds.




Figure 40. XICs of Buprenorphine at Various Concentrations in Urine (Diluted 10-fold, 10 µL injection)




Figure 41. XICs of Fentanyl at Various Concentrations in Urine (Diluted 10-fold, 10 µL injection)




Figure 42. XICs of Amitriptyline at Various Concentrations in Urine (Diluted 10-fold, 10 µL injection)

In the SWATH® Acquisition, MS/MS information is always available and so we can confirm the presence of the compound through MS/MS library matching (Figures 43 and 44) at the same time as determining how much of the compound is present.




Figure 43. Confident Identification of JWH-073 3-Butyl and Nortapentadol from SWATH® Acquisition Through Library Searching




Figure 44. Confident Identification of Amitriptyline and EDDP from SWATH® Acquisition Through Library Searching; Showing LC Separation Between Isomers was Still Achievable with this Fast Method

At the cutoff concentrations, library matching worked well with 80% of compounds yielding greater than 70% hit score (Figure 45).

Sample Name	Component Na...	Actual Conc...	Height	Ret. Time	Precursor Mass	Mass Error...	RT Conf...	Isotope Conf...	Library Conf...	Found At Mass	Mass Error...	Ret. Tim...	Library Hit	Library Score	Comb. Score	Isotope Ratio Diff.
EPDCCurve_UCa2	Alpha-Hydroxy...	50.00	39322	1.31	325.0811		✓	✓	✓	325.0811	0.0	0.00	Alpha-Hydroxy...	100.0	95.1	18.6
EPDCCurve_UCa2	Alpha-Hydroxy...	50.00	54742	1.32	342.0804		✓	✓	✓	342.0804	0.0	0.00	Alpha-Hydroxy...	100.0	99.3	1.8
EPDCCurve_UCa2	Buprenorphine	30.00	26422	1.24	468.1108		✓	✓	✓	468.1107	-0.3	0.01	Buprenorphine	100.0	97.9	2.9
EPDCCurve_UCa2	Desalofluazep...	50.00	39350	1.33	289.0538		✓	✓	✓	289.0539	0.2	0.01	Desalofluazep...	100.0	98.2	1.4
EPDCCurve_UCa2	EDDP	100.00	768379	1.20	278.1903		✓	✓	✓	278.1906	0.8	0.00	EDDP	100.0	96.1	1.3
EPDCCurve_UCa2	Fentanyl	2.00	10599	1.17	337.2274		✓	✓	✓	337.2276	0.5	0.01	Fentanyl	100.0	97.3	0.7
EPDCCurve_UCa2	Hydrocodone	50.00	61526	0.92	300.1594		✓	✓	✓	300.1597	0.9	0.01	Hydrocodone	100.0	94.9	3.3
EPDCCurve_UCa2	JWH-018 4-OH...	10.00	29593	1.39	358.1802		✓	✓	✓	358.1801	-0.1	0.01	JWH-018 Hydro...	100.0	98.0	4.6
EPDCCurve_UCa2	JWH-019 6-OH A...	10.00	44558	1.40	372.1958		✓	✓	✓	372.1959	0.3	0.01	JWH-019 N-Hyd...	100.0	96.7	4.4
EPDCCurve_UCa2	JWH-073 3-OH...	10.00	31459	1.38	344.1645		✓	✓	✓	344.1648	0.9	0.01	JWH-073 N-Hyd...	100.0	94.3	5.1
EPDCCurve_UCa2	JWH-122 5-OH...	10.00	44558	1.40	372.1958		✓	✓	✓	372.1959	0.3	0.01	JWH-019 N-Hyd...	100.0	96.7	4.4
EPDCCurve_UCa2	Martingapham...	100.00	13875	0.92	150.1277		✓	✓	✓	150.1277	-0.4	0.00	Martingapham...	100.0	98.0	0.8
EPDCCurve_UCa2	Methamphetamine	30.00	193407	1.07	234.1489		✓	✓	✓	234.1490	0.8	0.00	Methamphetamine	100.0	96.4	0.7
EPDCCurve_UCa2	Norbuprenorph...	30.00	27123	1.17	414.2639		✓	✓	✓	414.2641	0.5	0.00	Norbuprenorph...	100.0	97.5	2.3
EPDCCurve_UCa2	Nordiazepam	50.00	44694	1.35	271.0633		✓	✓	✓	271.0635	0.8	0.00	Nordiazepam	100.0	95.5	3.7
EPDCCurve_UCa2	Noroxycodone	50.00	23823	0.90	302.1387		✓	✓	✓	302.1386	-0.2	0.01	Noroxycodone	100.0	97.9	4.3
EPDCCurve_UCa2	Tapentadol	50.00	203639	1.09	222.1852		✓	✓	✓	222.1854	0.7	0.00	Tapentadol	100.0	96.0	1.1
EPDCCurve_UCa2	Tamoxifen	50.00	58735	1.34	301.0738		✓	✓	✓	301.0738	-0.2	0.01	Tamoxifen	100.0	98.1	3.2
EPDCCurve_UCa2	Tamozol	50.00	169458	1.06	264.1958		✓	✓	✓	264.1961	1.0	0.00	Tamozol	100.0	95.9	0.6
EPDCCurve_UCa2	Zolpidem	50.00	277586	1.33	308.1757		✓	✓	✓	308.1760	0.8	0.00	Zolpidem	100.0	96.1	1.9
EPDCCurve_UCa2	7-Aminoocozain...	50.00	76113	1.12	286.0742		✓	✓	✓	286.0744	0.9	0.01	7-Aminoocozain...	99.8	95.0	3.2
EPDCCurve_UCa2	Alprazolam	50.00	123684	1.32	309.0902		✓	✓	✓	309.0904	0.8	0.00	Alprazolam	99.8	95.7	4.0
EPDCCurve_UCa2	Diazepam	50.00	110405	1.36	285.0789		✓	✓	✓	285.0791	0.7	0.01	Diazepam	99.7	95.5	3.4
EPDCCurve_UCa2	Naloxone	50.00	56339	0.91	342.1700		✓	✓	✓	342.1698	-0.4	0.00	Naloxone	99.7	97.3	3.0
EPDCCurve_UCa2	Oxycodone	50.00	6763	0.41	302.1387		✓	✓	✓	302.1389	0.8	0.02	Oxycodone	99.7	94.2	3.5
EPDCCurve_UCa2	Amphetamine	100.00	3665	0.88	116.1121		✓	✓	✓	116.1120	-0.3	0.00	Amphetamine	99.6	97.1	4.1

Figure 45. Library Searching and Identification of Compounds in the 2.0 Minute Method at Cutoff Concentration Levels

MRM^{HR}

MRM^{HR} is a purely targeted data MS/MS acquisition and can be unscheduled or scheduled. The only non-targeted and therefore retrospective capability is through the TOF-MS experiment which is performed at the beginning of every scan. The power of the workflow however, is its selectivity capabilities through the accurate mass of unique fragment ions for quantification purposes. This is demonstrated in Figure 46 where MRM^{HR} is compared to the MRM analysis, extracted at nominal mass, and the extraction of the accurate mass of the precursor ion from a TOF-MS experiment. The compound is not able to be distinguished from the high background and interferences of the nominal mass experiment and not even by the extraction of the accurate mass of the precursor ion from the full scan TOF-MS experiment. It is not until we extract out two unique accurate mass fragment ions from the IMRM^{HR} experiment that we achieve the selectivity required to detect this compound by removal of the background and interferences and increase the S/N; improving the quantification capabilities. Another example of this selectivity gain over the accurate mass of the precursor ion is demonstrated in Figure 47 where a visible improvement in S/N is gained for the analysis of buprenorphine by the MRM^{HR} approach.




Figure 46. Increased Selectivity with MRM^{HR}; Avoiding False Negatives. Example given is a Feed Sample Tested Positive for NP Semicarbazide




Figure 39. Scheduled MRM^{HR} Selectivity Compared to TOF-MS; Buprenorphine (5ng/mL in urine, 10 fold dilution, 10 µL injection)

Quantification performance of the MRM^{HR} is demonstrated in Figure 40 for the 8.0 minute LC-MS/MS method.

Scheduled MRM^{HR}, 372.2 → 169.0644 ± 0.0100 m/z




Figure 40. MRM^{HR} Quantification Results for JWH-122 5-OH Pentyl in urine (Urine was diluted 10-fold, 10 µL injection)

Negative Mode Performance

Figures 41 and 42 show a couple of examples of negative mode performance of the SCIEX X500R QTOF System.



Figure 41. Negative Mode Performance of SCIEX X500R QTOF System for Analysis of Amobarbital



Figure 42. Negative Mode Performance of SCIEX X500R QTOF System for Analysis of THC-COOH

Conclusion

The arrival of the next generation QTOF, with the launch of the SCIEX X500R QTOF System and SCIEX OS Software, brings the powerful performance capabilities of the high resolution accurate mass technology to the routine identification and quantification forensic workflows.

- Hardware
 - SCIEX ExionLCTM Systems
 - Fully controlled by SCIEX OS software
 - Improved software integration for better stability
 - SCIEX X500R QTOF System
 - N-geometry design (same effective flight path length for ions and therefore resolution than V-geometry, but in a smaller overall footprint)
 - Heated TOF path for mass accuracy stability

- Minimized footprint, engineered for simplicity and service accessibility
 - Software
 - SCIEX OS Software
 - Intuitive and logical single software platform for LC control, MS control, data processing and reporting.
 - New user interface
 - Simultaneous identification and quantitation
- runtime with MS/MS information always being available with this MS/MS^{ALL} approach.

We have described the screening and quantification workflows of the SCIEX X500R QTOF System. Each workflow is straightforward to setup in the newly designed SCIEX OS Software and depending on the end users requirements we have demonstrated in this technical note the strengths of each workflow. Each provides TOF-MS and TOF-MS/MS analysis, both data being crucial in confidently identifying and quantifying forensic compounds.

- TOF-MS
- TOF-MS/MS
 - IDA
 - Non-targeted data acquisition
 - MS quantitation
 - Highest confidence screening with MS/MS information
 - MRM^{HR}
 - Targeted data acquisition for quantitation purpose
 - Can be performed unscheduled or scheduled
 - SWATH[®] Acquisition (with variable windows)
 - Non-targeted data acquisition
 - MS/MS for everything all the time
 - Screening and quantitation (MS/MS)
 - Library Searching and Ion Ratio

We evaluated different LC runtime methods. The longer method aided eluting all analytes throughout the entire gradient as evenly as possible in order to maximize triggering IDA MS/MS for all components and reduce the MRM^{HR} concurrency for quality of data (*Scheduled* MRM^{HR}). The library searching worked well for the SWATH[®] Acquisition in the 2.0 minute LC

For Research Use Only. Not for use in diagnostic procedures.

© 2016 AB SCIEX. AB Sciex is doing business as Sciex. The trademarks mentioned herein are the property of AB Sciex Pte. Ltd. or their respective owners. AB SCIEX™ is being used under license.

Document number: RUO-MKT-02-3786-A



Headquarters
500 Old Connecticut Path | Framingham, MA 01701 USA
Phone 508-383-7700
www.sciex.com

International Sales
For our office locations please call the division
headquarters or refer to our website at
www.sciex.com/offices