

Using the iDQuant™ Standards Kit for Pesticide Analysis to Analyze Residues in Fruits and Vegetable Samples

André Schreiber¹, Marjorie Minkoff², and Brent Lefebvre¹

¹ AB SCIEX, Concord, ON, Canada; ² AB SCIEX, Framingham, MS, USA

Overview

The iDQuant™ Standards Kit for Pesticide Analysis of certified reference material (ISP Guide 34, IOS/IEC 17025 and ISO 9001:2008) includes 204 pesticides. This kit has been designed to eliminate the need to source pesticides individually, and measure each one manually. This allows to setup, verify, and validate the performance of your pesticide screening method faster.

Here we present example data where we used the iDQuant™ Kit to screen for, quantify, and identify pesticides in extracts of fruits and vegetables using Liquid Chromatography tandem Mass Spectrometry (LC-MS/MS) with an AB SCIEX QTRAP® 5500 system.

Introduction

Many recent developments such as generic extraction procedures, (like QuEChERS), generic LC separation methods, and highly selective and sensitive MS/MS detectors have established LC-MS/MS as a standard technique for the analysis of pesticide residues in food. Hundreds of targeted analytes can be detected in a single analytical run. Software designed for ease-of-use like Cliquid® software in conjunction with iMethod™ applications made the fast adaptation of LC-MS/MS in new routine testing laboratories possible. For a while, data processing was the bottleneck for most laboratories, but the development of fast and automatic data processing and reporting tools sped up the delivery of analytical results.

A major hurdle for each laboratory is method setup, verification, and validation. The newly available iDQuant™ kit can be used for the following tasks:

- Tune or verify MRM transitions for best selectivity and sensitivity of MS/MS detection
- Measure retention times to quickly update acquisition methods using the *Scheduled MRM™* algorithm
- Verify the performance of LC methods
- Investigate recovery and reproducibility of sample preparation procedures
- Validate the performance of the complete method procedure



Description of the iDQuant™ Standards Kit for Pesticide Analysis

Content

The iDQuant™ Kit contains a total of 204 pesticides in 10 mixes at a concentration of 100 µg/mL and a certificate of analysis for each mixture. The pesticides are listed in Table 1.




Figure 1. Detection of 204 pesticides of the iDQuant™ Kit using a 4000 QTRAP® system and 2 MRM transitions per pesticide

Table 1. Pesticides in the iDQuant™ Standard Kit for Pesticide Analysis

Pesticide	CAS	Formula
Standard A1		
<i>3-Hydroxycarbofuran</i>	16655-82-6	C12H15NO4
<i>Acephate</i>	30560-19-1	C4H10NO3PS
<i>Aldicarb sulfone</i>	1646-88-4	C7H14N2O4S
<i>Aldicarb sulfoxide</i>	1646-87-3	C7H14N2O3S
<i>Aminocarb</i>	2032-59-9	C11H16N2O2
<i>Butocarboxim</i>	34681-10-2	C7H14N2O2 S
<i>Butoxycarboxim</i>	34681-23-7	C7H14N2O4S
<i>Carbendazim</i>	10605-21-7	C9H9N3O2
<i>Cymoxanil</i>	57966-95-7	C7H10N4O3
<i>Dicrotophos</i>	141-66-2	C8H16NO5P
<i>Dimethoate</i>	60-51-5	C5H12NO3PS2
<i>Dioxacarb</i>	6988-21-2	C11H13NO4
<i>Formetanate HCL</i>	23422-53-9	C11H15N3O2 . HCl
<i>Fuberidazole</i>	3878-19-1	C11H8N2O
<i>Methamidophos</i>	10265-92-6	C2H8NO2PS
<i>Methomyl</i>	16752-77-5	C5H10N2O2S
<i>Mevinphos</i>	7786-34-7	C7H13O6P
<i>Monocrotophos</i>	6923-22-4	C7H14NO5P
<i>Omethoate</i>	1113-02-6	C5H12NO4PS
<i>Oxamyl</i>	23135-22-0	C7H13N3O3S
<i>Propamocarb</i>	24579-73-5	C9H20N2O2
<i>Thiabendazole</i>	148-79-8	C10H7N3S
<i>Tricyclazole</i>	41814-78-2	C9H7N3S
<i>Vamidotion</i>	2275-23-2	C8H18NO4PS2
Standard B1		
<i>Aldicarb</i>	116-06-3	C7H14N2O2S
<i>Bendiocarb</i>	22781-23-3	C11H13NO4
<i>Carbaryl</i>	63-25-2	C12H11NO2
<i>Carbetamide</i>	16118-49-3	C12H16N2O3
<i>Carbofuran</i>	1563-66-2	C12H15NO3
<i>Carboxin</i>	5234-68-4	C12H13NO2S
<i>Clethodim</i>	99129-21-2	C17H26ClNO3S
<i>Desmedipham</i>	13684-56-5	C16H16N2O4
<i>Diethofencarb</i>	87130-20-9	C14H21NO4
<i>Ethiofencarb</i>	29973-13-5	C11H15NO2S

<i>Furalaxyl</i>	57646-30-7	C17H19NO4
<i>Iprovalicarb</i>	140923-17-7	C18H28N2O3
<i>Isoproc carb</i>	2631-40-5	C11H15NO2
<i>Metalaxyl</i>	57837-19-1	C15H21NO4
<i>Methiocarb</i>	2032-65-7	C11H15NO2S
<i>Mexacarbate</i>	315-18-4	C12H18N2O2
<i>Oxadixyl</i>	77732-09-3	C14H18N2O4
<i>Pirimicarb</i>	23103-98-2	C11H18N4O2
<i>Promecarb</i>	2631-37-0	C12H17NO2
<i>Propham</i>	122-42-9	C10H13NO2
<i>Propoxur</i>	114-26-1	C11H15NO3
<i>Pyracarbolid</i>	24691-76-7	C13H15NO2
<i>Thiofanox</i>	39196-18-4	C9H18N2O2S
<i>Thiophanate-methyl</i>	23564-05-8	C12H14N4O4S2
Standard C1		
<i>Alanycarb</i>	83130-01-2	C17H25N3O4S2
<i>Amitraz</i>	33089-61-1	C19H23N3
<i>Benalaxyl</i>	71626-11-4	C20H23NO3
<i>Benfuracarb</i>	82560-54-1	C20H30N2O5S
<i>Benzoximate</i>	29104-30-1	C18H18ClNO5
<i>Bifenazate</i>	149877-41-8	C17H20N2O3
<i>Cyazofamid</i>	120116-88-3	C13H13ClN4O2S
<i>Fenamidone</i>	161326-34-7	C17H17N3OS
<i>Fenazaquin</i>	120928-09-8	C20H22N2O
<i>Fenhexamid</i>	126833-17-8	C14H17Cl2NO2
<i>Fenoxycarb</i>	79127-80-3	C17H19NO4
<i>Flufenacet</i>	142459-58-3	C14H13F4N3O2S
<i>Furathiocarb</i>	65907-30-4	C18H26N2O5S
<i>Indoxacarb</i>	144171-61-9	C22H17ClF3N3O7
<i>Mefenacet</i>	73250-68-7	C16H14N2O2S
<i>Mepronil</i>	55814-41-0	C17H19NO2
<i>Piperonyl butoxide</i>	51-03-6	C19H30O5
<i>Quinoxifen</i>	124495-18-7	C15H8Cl2FNO
<i>Spiroxamine</i>	118134-30-8	C18H35NO2
<i>Thiobencarb</i>	28249-77-6	C12H16ClNOS
<i>Zoxamide</i>	156052-68-5	C14H16Cl3NO2

Standard D1

<i>Acetamiprid</i>	135410-20-7	C10H11ClN4
<i>Acibenzolar-S-methyl</i>	135158-54-2	C8H6N2OS2
<i>Bromucanazole</i>	116255-48-2	C13H12BrCl2N3O
<i>Clothianidin</i>	210880-92-5	C6H8ClN5O2S
<i>Cyproconazole</i>	113096-99-4	C15H18ClN3O
<i>Epoxiconazole</i>	135319-73-2	C17H13ClFN3O
<i>Etaconazole</i>	60207-93-4	C14H15Cl2N3O2
<i>Fenarimol</i>	60168-88-9	C17H12Cl2N2O
<i>Flutriafol</i>	76674-21-0	C16H13F2N3O
<i>Imazalil</i>	35554-44-0	C14H14Cl2N2O
<i>Imidacloprid</i>	138261-41-3	C9H10ClN5O2
<i>Metribuzin</i>	21087-64-9	C8H14N4OS
<i>Myclobutanil</i>	88671-89-0	C15H17ClN4
<i>Nitenpyram</i>	120738-89-8	C11H15ClN4O2
<i>Nuarimol</i>	63284-71-9	C17H12ClFN2O
<i>Paclobutrazol</i>	76738-62-0	C30H40Cl2N6O2
<i>Pyrimethanil</i>	53112-28-0	C12H13N3
<i>Thiacloprid</i>	111988-49-9	C10H9ClN4S
<i>Thiamethoxam</i>	153719-23-4	C8H10ClN5O3S
<i>Triadimenol</i>	55219-65-3	C14H18ClN3O2
<i>Triticonazole</i>	131983-72-7	C17H20ClN3O

Standard E1

<i>Abamectin</i>	71751-41-2	C48H72O14; C47H70O14
<i>Bitertanol</i>	55179-31-2	C20H23N3O2
<i>Bupirimate</i>	41483-43-6	C13H24N4O3S
<i>Butafenacil</i>	134605-64-4	C20H18ClF3N2O6
<i>Clofentezine</i>	74115-24-5	C14H8Cl2N4
<i>Cyprodinil</i>	121552-61-2	C14H15N3
<i>Diclobutrazol</i>	75736-33-3	C15H19Cl2N3O
<i>Difenoconazole</i>	119446-68-3	C19H17Cl2N3O3
<i>Diniconazole</i>	83657-24-3	C15H17Cl2N3O
<i>Ethofumesate</i>	26225-79-6	C13H18O5S
<i>Fenbuconazole</i>	114369-43-6	C19H17ClN4
<i>Flusilazole</i>	85509-19-9	C16H15F2N3Si
<i>Flutolanil</i>	66332-96-5	C17H16F3NO2
<i>Hexaconazole</i>	79983-71-4	C14H17Cl2N3O

<i>Mepanipyrim</i>	110235-47-7	C14H13N3
<i>Metconazole</i>	125116-23-6	C17H22ClN3O
<i>Methoxyfenozide</i>	161050-58-4	C22H28N2O3
<i>Penconazole</i>	66246-88-6	C13H15Cl2N3
<i>Prochloraz</i>	67747-09-5	C15H16Cl3N3O2
<i>Propiconazole</i>	60207-90-1	C15H17Cl2N3O2
<i>Rotenone</i>	83-79-4	C23H22O6
<i>Tebufenozide</i>	112410-23-8	C22H28N2O2
<i>Triflumizole</i>	68694-11-1	C15H15ClF3N3O

Standard F1

<i>Carfentrazone-ethyl</i>	128639-02-1	C15H14Cl2F3N3O3
<i>Doramectin</i>	117704-25-3	C50H74O14
<i>Emamectin-benzoate</i>	155569-91-8	C49H75NO13; C48H73NO13
<i>Eprinomectin</i>	123997-26-2	C50H75NO14; C49H73NO14
<i>Fluquinconazole</i>	136426-54-5	C16H8Cl2FN5O
<i>Hexythiazox</i>	78587-05-0	C17H21ClN2O2S
<i>Hydramethylnon</i>	67485-29-4	C25H24F6N4
<i>Ipconazole</i>	125225-28-7	C18H24ClN3O
<i>Ivermectin</i>	70288-86-7	C48H74O14
<i>Mesotrione</i>	104206-82-8	C14H13NO7S
<i>Moxidectin</i>	113507-06-5	C37H53NO8
<i>Propargite</i>	2312-35-8	C19H26O4 S
<i>Spinosad</i>	168316-95-8	C41H65NO10
<i>Spirodiclofen</i>	148477-71-8	C21H24Cl2O4
<i>Spiromesifen</i>	283594-90-1	C23H30O4
<i>Tebuconazole</i>	107534-96-3	C16H22ClN3O
<i>Tetraconazole</i>	112281-77-3	C13H11Cl2F4N3O

Standard G1

<i>Chlortoluron</i>	15545-48-9	C10H13ClN2O
<i>Cycluron</i>	2163-69-1	C11H22N2O
<i>Diuron</i>	330-54-1	C9H10Cl2N2O
<i>Fenuron</i>	101-42-8	C9H12N2O
<i>Fluometuron</i>	2164-17-2	C10H11F3N2O
<i>Forchlorfenuron</i>	68157-60-8	C12H10ClN3O
<i>Isoproturon</i>	34123-59-6	C12H18N2O

<i>Methabenzthiazuron</i>	18691-97-9	C10H11N3OS
<i>Methoprotryne</i>	841-06-5	C11H21N5OS
<i>Metobromuron</i>	3060-89-7	C9H11BrN2O2
<i>Monolinuron</i>	1746-81-2	C9H11ClN2O2
<i>Prometon</i>	1610-18-0	C10H19N5O
<i>Pymetrozine</i>	123312-89-0	C10H11N5O
<i>Secbumeton</i>	26259-45-0	C10H19N5O
<i>Simetryn</i>	1014-70-6	C8H15N5S
<i>Sulfentrazone</i>	122836-35-5	C11H10Cl2F2N4O3S
<i>Tebuthiuron</i>	34014-18-1	C9H16N4OS
<i>Terbumeton</i>	33693-04-8	C10H19N5O
<i>Thidiazuron</i>	51707-55-2	C9H8N4OS
<i>Triadimefon</i>	43121-43-3	C14H16ClN3O2

Standard H1

<i>Ametryn</i>	834-12-8	C9H17N5S
<i>Azoxystrobin</i>	131860-33-8	C22H17N3O5
<i>Boscalid</i>	188425-85-6	C18H12Cl2N2O
<i>Chloroxuron</i>	1982-47-4	C15H15ClN2O2
<i>Diflubenzuron</i>	35367-38-5	C14H9ClF2N2O2
<i>Dimethomorph</i>	110488-70-5	C21H22ClNO4
<i>Dimoxystrobin</i>	149961-52-4	C19H22N2O3
<i>Famoxadon</i>	131807-57-3	C22H18N2O4
<i>Fipronil</i>	120068-37-3	C12H4Cl2F6N4OS
<i>Fludioxonil</i>	131341-86-1	C12H6F2N2O2
<i>Hexaflumuron</i>	86479-06-3	C16H8Cl2F6N2O3
<i>Kresoxim-methyl</i>	143390-89-0	C18H19NO4
<i>Linuron</i>	330-55-2	C9H10Cl2N2O2
<i>Neburon</i>	555-37-3	C12H16Cl2N2O
<i>Phenmedipham</i>	13684-63-4	C16H16N2O4
<i>Picoxystrobin</i>	117428-22-5	C18H16F3NO4
<i>Prometryne</i>	7287-19-6	C10H19N5S
<i>Siduron</i>	1982-49-6	C14H20N2O
<i>Terbutryne</i>	886-50-0	C10H19N5S
<i>Triflumuron</i>	64628-44-0	C15H10ClF3N2O3

Standard J1

<i>Buprofezin</i>	69327-76-0	C16H23N3OS
<i>Chlorfluazuron</i>	71422-67-8	C20H9Cl3F5N3O3
<i>Ethiprole</i>	181587-01-9	C13H9Cl2F3N4OS
<i>Etoazole</i>	153233-91-1	C21H23F2NO2
<i>Fenpropimorph</i>	67564-91-4	C20H33NO
<i>Fenpyroximate</i>	111812-58-9	C24H27N3O4
<i>Fluazinam</i>	79622-59-6	C13H4Cl2F6N4O4
<i>Flufenoxuron</i>	101463-69-8	C21H11ClF6N2O3
<i>Fluoxastrobin</i>	361377-29-9	C21H16ClFN4O5
<i>Lufenuron</i>	103055-07-8	C17H8Cl2F8N2O3
<i>Novaluron</i>	116714-46-6	C17H9ClF8N2O4
<i>Prothioconazole</i>	178928-70-6	C14H15Cl2N3OS
<i>Pyraclostrobin</i>	175013-18-0	C19H18ClN3O4
<i>Pyridaben</i>	96489-71-3	C19H25ClN2OS
<i>Pyriproxyfen</i>	95737-68-1	C20H19NO3
<i>Tebufenpyrad</i>	119168-77-3	C18H24ClN3O
<i>Teflubenzuron</i>	83121-18-0	C14H6Cl2F4N2O2
<i>Trifloxystrobin</i>	141517-21-7	C20H19F3N2O4

Standard K1

<i>Chlorantraniliprole</i>	500008-45-7	C18H14BrCl2N5O2
<i>Cyromazine</i>	66215-27-8	C6H10N6
<i>Dinotefuran</i>	165252-70-0	C7H14N4O3
<i>Ethirimol</i>	23947-60-6	C11H19N3O
<i>Fenobucarb</i>	3766-81-2	C12H17NO2
<i>Fonicamid</i>	158062-67-0	C9H6F3N3O
<i>Flubendimide</i>	272451-65-7	C23H22F7IN2O4S
<i>Halofenozide</i>	112226-61-6	C18H19ClN2O2
<i>Isocarbophos</i>	24353-61-5	C11H16NO4PS
<i>Mandipropamid</i>	374726-62-2	C23H22ClNO4
<i>Metaflumizone</i>	139968-49-3	C24H16F6N4O2
<i>Pencycuron</i>	66063-05-6	C19H21ClN2O
<i>Spinetoram</i>	187166-40-1	C42H69NO10
<i>Spirotetramat</i>	203313-25-1	C21H27NO5
<i>Temephos</i>	3383-96-8	C16H20O6P2S3

Storage and Stability

The pesticides of the iDQuant™ Kit are stable in the sealed ampoules for one year at -20°C. An example of stability test results is shown in Figure 2.




Figure 2. Results of the stability testing indicating stability of all pesticides over a period of one year at -20°C

Stock solutions of the 204 pesticides at 1 µg/kg can be simply prepared by mixing 10 µL of each mix and diluting with 900 µL of acetonitrile. It is recommended to store aliquots of this stock solution in amber glass vials at -20°C, for future use.

Calibration standards can be prepared by serial dilution of the prepared stock solutions using water or the aqueous mobile phase of the LC system. It is recommended to prepare calibration standards daily to avoid degradation of pesticides.

LC-MS/MS Method Setup

The iMethod™ application for pesticide screening contains a ready-to-use method to screen for all 204 pesticides of the iDQuant™ kit.

The kit can also be used to fine tune MRM conditions or to adjust retention times of existing methods to update the *Scheduled* MRM™ algorithm for best performance (Figure 3).

Application of the iDQuant™ Standards Kit for Pesticide Analysis

Experimental

- The iDQuant™ standards kit for pesticide analysis was used for method setup and preparation of calibration standards.




Figure 3. The injection of each mix and processing in MultiQuant™ software can be used to quickly update retention times in the *Scheduled* MRM™ method for best performance

- Fruit and vegetable samples were extracted using a QuEChERS procedure and diluted 10x with water to optimize chromatographic peak shape and minimize possible matrix effects and interferences.¹⁻²
- LC separation was achieved on a Shimadzu UFLC_{XR} system with a Restek Ultra Aqueous C18 3 µm (100x2.1 mm) column and a 15 min gradient of water and methanol with ammonium formate buffer at a flow rate of 0.5 mL/min. The injection volume was set to 10 µL.
- The AB SCIEX QTRAP® 5500 system was operated with Turbo V™ source and Electrospray Ionization (ESI) probe. 436 MRM transitions were monitored using the *Scheduled* MRM™ algorithm with a detection window of 120 s and a target scan time of 0.7 s.
- MultiQuant™ 2.0.2 software with the 'Multicomponent' query was used for quantitative data processing (Figure 4).

Table 2. Pesticides quantified in fruit and vegetable samples above 10 µg/kg and positively identified using their retention time and MRM ratio

Sample	Pesticide	Concentration (µg/kg)	Retention time (min)	Expected retention time (min)	MRM ratio	Expected MRM ratio
<i>Tomato (grape)</i>	Imidacloprid	16.1	4.44	4.50	1.02	0.98
	Methoxyfenozide	21.0	7.42	7.50	0.25	0.23
<i>Tomato (vine)</i>	Boscalid	33.8	7.55	7.60	0.30	0.34
	Pyraclostrobin	11.4	8.47	8.50	0.76	0.77
	Pyridaben	18.6	9.79	9.80	0.93	0.91
	Tebuconazole	42.4	7.90	7.90	0.40	0.43
<i>Raspberry</i>	Azoxystrobin	76.4	7.34	7.40	0.18	0.18
	Boscalid	42.0	7.52	7.60	0.31	0.34
	Myclobutanil	46.7	7.75	7.80	0.56	0.65
	Pyraclostrobin	10.4	8.44	8.50	0.80	0.77
<i>Strawberry</i>	Boscalid	23.0	7.54	7.60	0.30	0.34
	Cyprodinil	203	8.67	8.70	0.82	0.93
	Fenhexamid	141	7.59	7.60	0.94	0.87
	Fludioxinil	127	7.47	7.50	0.25	0.25
	Novaluron	27.8	8.82	8.80	0.70	0.84
	Pyrimethanil	24.2	7.81	7.90	0.21	0.22
<i>Banana</i>	Imazalil	15.0	8.58	8.60	0.16	0.19
<i>Lime</i>	Imazalil	230	8.60	8.60	0.17	0.19
	Thiabendazole	41.7	6.38	6.40	0.89	0.89
<i>Grapefruit</i>	Imazalil	82.2	8.60	8.60	0.24	0.19
	Thiabendazole	299	6.37	6.40	0.90	0.89
<i>Red pepper</i>	Azoxystrobin	12.4	7.37	7.40	0.16	0.18
	Clothianidin	18.4	4.15	4.20	1.05	1.14
	Imidacloprid	83.0	4.46	4.50	0.97	0.98
	Tebuconazole	41.1	7.90	7.90	0.36	0.43
	Thiamethoxam	64.0	3.61	3.60	0.84	0.96
<i>Grapes not washed</i>	Boscalid	15.6	7.55	7.60	0.34	0.34
	Dimethomorph	12.1	7.82	7.80	0.59	0.65
	Famoxadone	650	8.24	8.20	0.47	0.52
	Tebuconazole	130	8.32	8.4	0.14	0.11
<i>Grapes washed</i>	Boscalid	11.7	7.54	7.60	0.37	0.34
	Dimethomorph*	5.91	7.81	7.80	0.63	0.65
	Famoxadone	367	8.23	8.20	0.50	0.52
	Tebuconazole	60.3	8.31	8.40	0.13	0.11

* Dimethomorph in the grapes sample after washing was quantified and identified below 10 µg/kg




Figure 4. Result display after applying the 'Multicomponent' query in MultiQuant™ software, only pesticides above 10 µg/kg and with a positive identification based on their MRM ratio are displayed

Results

All identified and quantified pesticides are listed in Table 2. The 'Multicomponent' query was used in MultiQuant™ software to automatically report concentrations above the threshold of 10 µg/kg and to compare MRM ratios for compound identification. The maximum permitted tolerances for MRM intensities of the SANCO guideline 'Method validation and quality control procedures for pesticide residue analysis in food and feed' were used for identification.³

As part of the samples set grapes were analyzed with and without washing before extraction. Half of the grapes were simply rinsed under running tap water for approximately 30 seconds.

For Research Use Only. Not for use in diagnostic procedures.

© 2011 AB SCIEX. The trademarks mentioned herein are the property of AB Sciex Pte. Ltd. or their respective owners. AB SCIEX™ is being used under license.

Publication number: 3370211-01

The results in Table 2 show that this washing method is not sufficient to completely remove fungicide before consumption.

Summary

The iQuant™ Standards Kit for Pesticide Analysis was used to successfully setup an LC-MS/MS method to identify and quantify pesticides in fruit and vegetable samples. This was accomplished in less time than it was previously possible when a standard kit was either unavailable or had to be created from individual standards.

References

- 1 M. Anastassiades, et al.: J. AOAC Int. 86 (2003) 412-431
- 2 J. Wong et al.: J. Agric. Food Chem. 58 (2010) 5897-5903
- 3 SANCO/10684/2009