

Simultaneous detection and quantification of 15 drugs of abuse in whole blood by online solid phase extraction and LC-MS/MS

3200 QTRAP[®] LC/MS/MS system

H. Grison-Hernando,¹ C. Renaud,¹ A. Morla,² G. Deslandes,¹ A. Pineau,¹ E. Dailly,¹ P. Jolliet,¹ C. Monteil-Ganiere.¹

¹Laboratoire de Pharmacologie-Toxicologie, CHU de Nantes, France; ²AB Sciex, Les Ulis, France.

Introduction

A new method for quantification of drugs of abuse in whole blood was developed with a simple sample pre-treatment, online solid-phase extraction and liquid chromatography-tandem mass spectrometry (LC-MS/MS). Fifteen drugs of abuse and metabolites were measured in a single chromatographic run. These included amphetamine, metamphetamine, MDMA, MDA, MDEA, MBDB, mephedrone, 6-MAM, morphine, codeine, dihydrocodeine, ethylmorphine, cocaine, BEG and cocaethylene).

Materials and Methods

Sample Preparation

For quantitative determination, 14 deuterated analogues were used as internal standards. After protein precipitation of 250 μ L whole blood with ZnSO₄-methanol spiked with a mixture of internal standards deuterated, samples were mixed and centrifuged; the supernatant was evaporated to dryness, and reconstituted with mobile phase, before injection in the chromatographic system as described below.

- Added 250 μ L whole blood to a microcentrifuge tube
- Added 500 μ L ZnSO₄·7H₂O 0.2 M - MeOH (1:4) + 14 internal standards
- Vortexed followed by centrifugation
- Removed the supernatant and placed in HPLC vial and dried under nitrogen
- Added 100 μ L of mobile phase A

HPLC Conditions

Analytes were firstly loaded on the extraction column (Strata-X, Phenomenex) with a 2 mM ammonium formate buffer, diverting salts and unwanted components to waste. A switching valve was then triggered to back-flush with mobile phase and elute analytes from the extraction column and diverting the flow onto the analytical column (Kinetex PFP, Phenomenex) kept at 60°C. The compounds were well separated in a total run time of 15 minutes. A gradient was performed with the aqueous phase A (2 mM ammonium formate and 0.2% formic acid), and the organic phase B (2 mM ammonium formate with methanol-acetonitrile (70-30) and 0.2% formic acid) (figure 1).




Figure 1. Chromatographic gradient profile

MS/MS Conditions

The mass spectrometer used was an AB SCIEX 3200 QTRAP[®] LC/MS/MS System. The method was developed using the *scheduled* MRM[™] algorithm, with two transitions per compound. All analytes were detected in positive ionization mode (figure 2).

Analyte Internal Standard	Q1 (m/z)	Q3 (m/z)	DP (V)	CE (V)	Analyte Internal Standard	Q1 (m/z)	Q3 (m/z)	DP (V)	CE (V)
<i>Amphetamine</i>	136.1	91	24	25	<i>6-MAM</i>	328.2	165.1	56	45
		65		50			211.2		31
<i>D5-Amphetamine</i>	141	93.1	40	20	<i>D3-6-MAM</i>	331.1	165.1	56	53
<i>Metamphetamine</i>	150.1	91	40	27	<i>Morphine</i>	286.1	152.2	56	73
		65		50			165.2		47
<i>D5-Metamphetamine</i>	155	92.1	40	20	<i>D3-Morphine</i>	289	165.2	40	50
<i>MDMA</i>	194.1	163.1	31	17	<i>Codeine</i>	300.2	165.2	56	61
		105.1		34			215.2		34
<i>D5-MDMA</i>	199.1	165.1	31	17	<i>D3-Codeine</i>	303	215.2	40	35
<i>MDA</i>	180.1	105.1	21	27	<i>Dihydrocodeine</i>	302.2	199.2	43	46
		135		23			171.2		57
<i>D5-MDA</i>	185.1	110.1	88	30	<i>D6-Dihydrocodeine</i>	308.1	202.2	61	43
<i>MDEA</i>	208.1	105.1	26	33	<i>Ethylmorphine</i>	314.1	152.2	66	89
		163.2		17			115.1		101
<i>D5-MDEA</i>	213.1	163.2	76	17	<i>Cocaine</i>	304.2	182.2	66	89
<i>MBDB</i>	208	77	26	57			105.1		101
		51		79	<i>D3-Cocaine</i>	307	185.2	40	20
<i>D5-MBDB</i>	213.1	136.1	40	20					
<i>Mephedrone</i>	178.1	160.2	26	15	105.1	39			
		144.2		35	<i>D3-Benzoylcegonine</i>	293	171.2	40	20
<i>D3-Mephedrone</i>	181.1	148.1	26	23					
					82.1	43			
					<i>D3-Cocaethylene</i>	321.1	199.3	36	25

Figure 2. MRM transitions and compound- dependent parameters

Results

The dynamic range of the assay for 6-MAM and Cocaine was demonstrated to be from 0 to 50 ng/mL with linearity up to 250 ng/mL. All other compounds showed a dynamic range of 0-200 ng/mL, linear up to 1000 ng/mL. To assess precision and repeatability, 5 replicate injections of 3 QC levels were performed. For reproducibility, the extraction was performed on 5 separate days. For all 3 QC levels the coefficient of variations (CVs) were less than 15% and accuracies of between 85 and 115% for all experiments. Limit of quantitation for 6-MAM and Cocaine were both 1.25 ng/mL and all other compounds was 5 ng/mL.

Recoveries measured with addition before and after precipitation were all greater than 75%. Compounds after precipitation were shown to be stable for 72 hours at 10 °C and cross contamination was shown to be less than 0.6%. Matrix effect was evaluated using six different whole bloods at 2 concentration levels with a result of coefficient of variation of less than 15%. A good correlation was demonstrated when the results of the developed LC/MS/MS were compared to the results of a GC/MS method that consisted of a liquid/liquid extraction followed by derivitization.




Figure 3. Example Calibrator Chromatogram

Conclusion

An LC/MS/MS method has been developed that enables the simultaneous analysis of amphetamines, cocaine, opiates and mephedrone from whole blood. The method has a simple and fast sample preparation with a short chromatographic separation. It provides an alternative approach of the lengthy liquid/liquid extraction and derivitization procedure required for GC/MS analysis, using low volume of sample and with no matrix effects.

For Research Use Only. Not for use in diagnostic procedures.

© 2014 AB SCIEX. The trademarks mentioned herein are the property of AB Sciex Pte. Ltd. or their respective owners. AB SCIEX™ is being used under license.

RUO-MKT-02-1225-A