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Comparation the N-Glycan of different source's recombinant human erythropoietin by using capillary electrophoresis with laser-induced fluorescence detector

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¹ CE & Biopharma, Sciex, China

² Division of Recombinant Biological Products, National Institute of Food and Drug Control, Beijing, China

Key Words: rhuEPO, Glycosylation, CE-LIF, Fast

Abstract: Recombinant human erythropoietin (rhuEPO) has been extensively used as a pharmaceutical product for treating anemia in clinic. The rhuEPO is a typical glycoprotein, and the glycosylation is up to 40%. So, glycosylation of rhuEPO is of great importance and could affects the biological activity, immunogenicity and pharmacokinetics. In this study, by using effective glycan capture by magnetic microparticles, we can realize the glycan releasing and labeling in 60 mins, and analysis in conjunction with a rapid high-resolution capillary electrophoresis analysis method for the prepared samples using laser induced fluorescence detection (CE-LIF) in 5 mins. Also, the GU software can identify the glycan structure automatically. As far as we know, this is the first time to compare different company's glycan of rhuEPO qualitatively and quantitatively by using CE-LIF. The proposed method is fast and has good reproducibility, which allowed the inter-day RSD of migration time and peak area were around 0.2% and 2.6% (n=5), respectively. Hence, this method has great significance of studying the glycan glycosylation of rhuEPO.

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- [1] European Pharmacopoeia Erythropoietin concentrated solution Elaboration of a European Pharmacopoeia[S]. 1528-1532.
- [2] Li F, Progress in research on glycosylation analysis of erythropoietin, [J]. Chin J Biologicals, 2018, 31, (3): 328.
- [3] Fast Glycan Labeling And Analysis Kit Application Guide, Sciex, B94499 PTO .

For Research Use Only. Not for use in Diagnostics Procedures.

AB Sciex is operating as SCIEX.

© 2019. AB Sciex. The trademarks mentioned herein are the property of AB Sciex Pte. Ltd. or their respective owners. AB SCIEX™ is being used under license.

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