
---

# MultiQuant™ MD 3.0.3

Vodič za korisnički interfejs

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

Slika 1-2 Meni „File“

File	Edit	Process	Audit Trail	Window	Help
New Results Table...					Ctrl+N
New Quantitation Method...					
Open Results Table...					Ctrl+O
Open Quantitation Method...					
Save					Ctrl+S
Save As...					
Recent Results Tables					▶
Recent Quantitation Methods					▶
Export					▶
Transfer to LIMS...					Ctrl+L
Export and Save Results Table...					
Create Report and Save Results Table...					Ctrl+R
Exit					


Engleski jezik	Prevod
New Results Table	Nova tabela sa rezultatima
New Quantitation Method	Novi metod kvantifikacije
Open Results Table	Otvori tabelu sa rezultatima
Open Quantitation Method	Otvori metod kvantifikacije
Save	Save
Save As	Sačuvaj kao
Recent Results Tables	Nedavne tabele sa rezultatima
Recent Quantitation Methods	Nedavni metodi kvantifikacije
Export	Opcija „Export“
Transfer to LIMS	Prenesi u LIMS
Export and Save Results Table	Dijalog „Export and Save Results Table“
Create Report and Save Results Table	Kreiraj izveštaj i sačuvaj tabelu sa rezultatima
Exit	Izlaz

**Slika 1-3 Meni „Edit“**


Engleski jezik	Prevod
Clear	Obriši
Copy	Kopiraj
Paste	Zalepi
Copy Entire Table	Kopiraj celu tabelu
Fill Down	Popuni nadole
Select All Rows	Izaberi sve redove
Modify Results Table Method	Menjanje metoda tabele sa rezultatima
Project Integration Defaults	Podrazumevane vrednosti integracije projekta
Project Units & Calibration Defaults	Podrazumevane vrednosti jedinica i kalibracije projekta
Project Secure Export Settings	Podešavanja bezbednog izvoza projekta
Enable Project Modified Peak Warning	Omogući upozorenje na izmenjenu vršnu vrednost projekta
Cache Chromatograms for Faster Peak Review	Keširaj hromatograme za brži pregled vršnih vrednosti
Cache all Chromatograms Now	Keširaj odmah sve hromatograme

Slika 1-4 Meni „Process“


Engleski jezik	Prevod
Add Samples	Dodaj uzorke
Remove Selected Samples	Ukloni izabrane uzorke
Show Only Outliers	Prikaži samo odstupanja
Go to Next Outlier	Pređi na sledeće odstupanje
Export Calibration and Save Results Table	Izvezi kalibraciju i sačuvaj tabelu sa rezultatima
Import External Calibration	Uvezi eksternu kalibraciju
Remove External Calibration	Ukloni eksternu kalibraciju

Slika 1-5 Meni „Audit Trail“

Engleski jezik	Prevod
Audit Trail Viewer	Prikazivač traga revizije
Audit Map Manager	Upravljač mapom revizije
View Session Audit Map	Prikaži mapu revizije sesije

**Slika 1-6 Meni „Help“**


Engleski jezik	Prevod
Install License	Instaliraj licencu
Verify Installation	Proveri instalaciju
Software Reference Guide	Referentni vodič za softver
Software Release Notes	Software Release Notes
About	Više informacija

# Kreiranje nove tabele sa rezultatima „Results Table“

## 2

Slika 2-1 Dijalog „Create Results Table - Select Samples“


Engleski jezik	Prevod
Select the desired samples	Izaberite željene uzorke
Available	Dostupno
Selected	Izabrano
Back	Nazad
Next	Sledeće

## Kreiranje nove tabele sa rezultatima „Results Table“

Engleski jezik	Prevod
Finish	Završi
Cancel	Otkazi

Slika 2-2 Dijalog „Create Results Table - Select Method“


Engleski jezik	Prevod
Select an existing quantitation method or create a new method now.	Izaberite postojeći metod kvantifikacije ili kreirajte novi metod
Choose Existing Method	Izaberite postojeći metod
Method Name	Naziv metoda
Open	Otvori
Edit method	Uredi metod

## Kreiranje nove tabele sa rezultatima „Results Table“

Engleski jezik	Prevod
Create New Method	Kreiraj novi metod
Method Name	Naziv metoda
New	Novo
Back	Nazad
Next	Sledeće
Finish	Završi
Cancel	Otkazi


Slika 2-3 Dijalog „Save Quantitation Method As“

## Kreiranje nove tabele sa rezultatima „Results Table“

Engleski jezik	Prevod
Save in	Sačuvaj u
File name	Naziv datoteke
Save	Sačuvaj
Cancel	Otkazi

Slika 2-4 Dijalog „Create Results Table - Select Representative Sample“


Engleski jezik	Prevod
The recommended Representative Sample has been automatically selected. Change the selection only if required.	Preporučeni reprezentativni uzorak je automatski izabran. Promenite izbor samo ako je to potrebno.
Sample Name	Naziv uzorka
Type	Tip

## Kreiranje nove tabele sa rezultatima „Results Table“

Engleski jezik	Prevod
Back	Nazad
Next	Sledeće
Finish	Završi
Cancel	Otkazi

Slika 2-5 Dijalog „Create Results Table - Define Components“


Engleski jezik	Prevod
Select or verify the analyte and internal standard names and masses.	Izaberite ili potvrdite nazive i mase analita i internih standarda.
Experiment	Eksperiment
Row	Red

## Kreiranje nove tabele sa rezultatima „Results Table“

Engleski jezik	Prevod
IS (internal standard)	IS (interni standard)
Name	Naziv
Group	Grupa
IS Name	Naziv IS-a
Back	Nazad
Next	Sledeće
Finish	Završi
Cancel	Otkazi

**Slika 2-6 Kontekstualni meni tabele „Create Results Table - Define Components Groups“**


Engleski jezik	Prevod
Clear	Obriši
Copy	Kopiraj
Paste	Zalepi
Find Component by Name	Pronađi komponentu prema nazivu
Insert Row Above	Umetni red iznad
Delete Selected Rows	Obriši izabrane redove
Sum Multiple Ions	Saberi višestruke jone

## Kreiranje nove tabele sa rezultatima „Results Table“

Engleski jezik	Prevod
Groups	Grupe
Using Constant Group Size	Korišćenje konstantne veličine grupe
By filling Down Existing Groups	Popunjavanjem postojećih grupa
Using Q1 Masses	Korišćenje Q1 masa
Using Q3 Masses	Korišćenje Q3 masa
Using (Q1 - Q3) Mass Differences	Korišćenje razlika masa (Q1 - Q3)
Add Group to Start of Component Name	Dodaj grupu na početak naziva komponente
Remove Group from Start of Component Name	Ukloni grupu sa početka naziva komponente
Append Summed Ions for Groups	Dodaj sabrane jone za grupe

**Slika 2-7 Kontekstualni meni tabele „Create Results Table - Define Components Internal Standards“**


Engleski jezik	Prevod
Clear	Obriši
Copy	Kopiraj
Paste	Zalepi
Find Component by Name	Pronađi komponentu prema nazivu
Insert Row Above	Umetni red iznad

## Kreiranje nove tabele sa rezultatima „Results Table“

Engleski jezik	Prevod
Delete Selected Rows	Obriši izabrane redove
Sum Multiple Ions	Saberi višestruke jone
Internal Standards	Interni standardi
Set IS for All Analytes	Postavi IS za sve analite
Set IS for Selected Analytes	Postavi IS za izabrane analite
Set Last Component of Group as IS	Postavi poslednju komponentu grupa kao IS
Set for All Groups as for Selected Group	Postavi za sve grupe kao za izabranu grupu

Slika 2-8 Dijalog „Create Results Table - Define Integration“


## Kreiranje nove tabele sa rezultatima „Results Table“

Engleski jezik	Prevod
Scroll through the various components, setting the integration parameters as needed.	Pomerajte prikaz kroz razne komponente i postavite parametre integracije po potrebi.
Gaussian Smooth Width: points	Širina Gausovog ravnjanja: tačke
Expected RT: min	Očekivano vreme zadržavanja: min
RT Half Window: sec	Poluprozor očekivanog vremena zadržavanja: sek
Updated Expected RT: No or Yes	Ažurirano očekivano vreme zadržavanja: ne ili da
Report Largest peak	Prijavi najveću vršnu vrednost
Min. Peak Width: points	Min. širina vršne vrednosti: tačke
Min Peak Height:	Min. visina vršne vrednosti:
Integration Parameters	Parametri integracije
Noise Percentage	Procenat šuma
Baseline Sub. Window: min	Prozor za skidanje osnovne linije: min
Peak Splitting: points	Deljenje vršne vrednosti: tačke
Regression	Regresija
Parameter: Area or Height	Parametar: površina ili visina
Fit: Linear, Linear Through Zero, Mean Response Factor, Quadratic, Power, Wagner, Hill	Okvir: linearni, linearni do nule, faktor srednje vrednosti odgovora, kvadratni, stepenovani, „Wagner“, „Hill“
Weighting: None, 1/x, 1x <sup>2</sup> , ln (x), 1/y, 1y <sup>2</sup> , ln (y)	Ponderisanje: nema, 1/x, 1x <sup>2</sup> , ln (x), 1/y, 1y <sup>2</sup> , ln (y)
Concentration Units	Jedinice koncentracije
Apply units to all Analytes	Primeni na sve analite
Back	Nazad
Next	Sledeće
Finish	Završi
Cancel	Otkazi

## Kreiranje nove tabele sa rezultatima „Results Table“


---

**Slika 2-9 Kontekstualni meni tabele „Create Results Table - Define Integration“**


Engleski jezik	Prevod
Find Component by Name	Pronađi komponentu prema nazivu
Highlight Components with Uncertain RT	Označi komponente sa neizvesnim vremenom zadržavanja
Home Graph Axis	Osa početnog grafikona
Overlay Other Components for Group	Preklopi ostale komponente za grupu
Update Retention Times	Dijalog „Update Retention Times“

**Slika 2-10 Dijalog „Find Component by Name“**


Engleski jezik	Prevod
Select component containing the following (case insensitive) text:	Izaberi komponente koje sadrže sledeći tekst (nema razlike između malih i velikih slova):
Select text:	Izaberi tekst:
OK	U redu
Cancel	Otkazi

Slika 2-11 Dijalog „Highlight Components“

Engleski jezik	Prevod
Highlight components with peak are less than _% of the total area.	Označi komponente čija je vršna vrednost manja od _% ukupne površine.
Only show this dialog again if the shift key is down	Prikaži dijalog ponovo samo kada je pritisnut taster Shift
OK	U redu
Cancel	Otkazi

Slika 2-12 Dijalog „Update Retention Times“

Engleski jezik	Prevod
Update retention times using the biggest peak in a window centered at the expected RT of width _ times the RT window	Ažuriraj vreme zadržavanja pomoću najveće vršne vrednosti u periodu centriranom na očekivano vreme zadržavanja za širinu _ puta period vremena zadržavanja
Only show this dialog again if the shift key is down	Prikaži dijalog ponovo samo kada je pritisnut taster Shift
OK	U redu
Cancel	Otkazi

Slika 2-13 Dijalog „Create Results Table - Outlier Settings“

MQ Create Results Table - Outlier Settings

Set criteria for flagging outliers.

☒ Accuracy for Standards

Max. Accuracy Tolerance for LLOQ (lowest Std):  %

Max. Accuracy Tolerance for Stds except LLOQ:  %

☐ Ion Ratio

☒ Calculated Concentration

Component	Lower Limit of Calculated Concentration	Upper Limit of Calculated Concentration
Flunitrazepam.IS		
Nordiazepam.IS		
Flunitrazepam		
Nordiazepam		

< Back   Next >   Finish   Cancel

Engleski jezik	Prevod
Set criteria for flagging outliers.	Postavi kriterijume za označavanje odstupanja
Accuracy for Standards	Preciznost za standarde
Max. Accuracy Tolerance for LLOQ (lowest Std)	Maks. tolerancija preciznosti za LLOQ (najniži standard)
Max. Accuracy Tolerance for Stds except LLOQ	Maks. tolerancija preciznosti za standarde osim LLOQ
Accuracy for QCs	Preciznost za kontrolu kvaliteta
Max. Accuracy Tolerance for QC	Maks. tolerancija preciznosti za kontrolu kvaliteta
Ion Ratio	Odnos jona

## Kreiranje nove tabele sa rezultatima „Results Table“

Engleski jezik	Prevod
Calculated Concentration	Izračunata koncentracija
Component	Komponenta
Lower Limit of Calculated Concentration	Donja granica izračunate koncentracije
Upper Limit of Calculated Concentration	Gornja granica izračunate koncentracije
Back	Nazad
Next	Sledeće
Finish	Završi
Cancel	Otkazi

Slika 2-14 Kontekstualni meni tabele „Create Results Table - Outlier Settings“


Component	Lower Limit of Calculated Concentration	Upper Limit of Calculated Concentration
► Flunitrazepam.IS		
Nordiazepam.IS		
Flunitrazepam		
Nordiazepam		
<div> <p>Apply to all analytes the Lower Limit of Calc. Concentration</p> <p>Apply to all analytes the Upper Limit of Calc. Concentration</p> </div>		

Engleski jezik	Prevod
Apply to all analytes the Lower Limit of Calc. Concentration	Primeni na sve analite za donju granicu izrač. koncentracije
Apply to all analytes the Upper Limit of Calc. Concentration	Primeni na sve analite za gornju granicu izrač. koncentracije

# Novi metod kvantifikacije


# 3

Slika 3-1 Kontekstualni meni „New Quantitation Method – Integration & Regression“


Engleski jezik	Prevod
Find Component by Name	Pronađi komponentu prema nazivu
Highlight components with Uncertain RT	Označi komponente sa neizvesnim vremenom zadržavanja
Home Graph Axis	Osa početnog grafikona
Overlay Other Components for Group	Preklopi ostale komponente za grupu
Update Retention Times	Ažuriraj vreme zadržavanja
Set New Typical Sample	Postavi novi tipičan uzorak

Slika 3-2 Dijalog „Select by Name“


Engleski jezik	Prevod
Select component containing the following (case insensitive) text:	Izaberi komponente koje sadrže sledeći tekst (nema razlike između malih i velikih slova):
Select text	Izaberi tekst
OK	U redu
Cancel	Otkazi

**Slika 3-3 Dijalog „Highlight Components“**

Engleski jezik	Prevod
Highlight components with peak area less than _% of the total area.	Označi komponente čija je oblast vršne vrednosti manja od _% ukupne oblasti.
Only show this dialog again if the shift key is down	Prikaži dijalog ponovo samo kada je pritisnut taster Shift
OK	U redu
Cancel	Otkazi

**Slika 3-4 Dijalog „Update Retention Times“**

## Novi metod kvantifikacije

---


Engleski jezik	Prevod
Update retention times using the biggest peak in a window centered at the expected RT of width _ times the RT window	Ažuriraj vreme zadržavanja pomoću najveće vršne vrednosti u periodu centriranom na očekivano vreme zadržavanja za širinu _ puta period vremena zadržavanja
Only show this dialog again if the shift key is down	Prikaži dijalog ponovo samo kada je pritisnut taster Shift
OK	U redu
Cancel	Otkazi

# Opcija „Export“


# 4

Slika 4-1 Dijalog „Export Metric“











Engleski jezik	Prevod
Metric: Area, IS Area, Corrected Area, IS Corrected Area, Area Ratio, Height, IS Height, Corrected Height, IS Corrected Height, Height Ratio, Area/Height, IS Area/Height, Corrected Area/Height, IS Corrected Area/Height, Region Height, IS Region Height, Quality, IS Quality, Retention Time, IS Retention Time,	Metrika: Oblast, IS oblast, Oblast ispravke, IS oblast ispravke, Odnos oblasti, Visina, IS visina, Ispravljena visina, IS ispravljena visina, Odnos visina, Oblast/Visina, IS oblast/visina, Ispravljena oblast/visina, IS ispravljena oblast/visina, Visina regiona, IS visina regiona, Kvalitet, IS kvalitet, Vreme zadržavanja, IS vreme zadržavanja
Only show this dialog again if the shift key is down	Prikaži dijalog ponovo samo kada je pritisnut taster Shift
OK	U redu
Cancel	Otkazi

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© 2017. AB Sciex