



## Advanced PA 800 Plus Troubleshooting Course (PN 149850)

Language: German or English

Training dates and Locations:

17-18 May 2016 SCIEX training facility Warrington, UK

6-7 September 2016 SCIEX training facility Darmstadt, Germany

### Day1 10AM-5.30PM

Sciex CE Plus family Systems Overview

Principles of Capillary Electrophoresis (presentation and interactive session)

General technical session

Selected Applications: FZCE SDS-MW/IgG Assay, cIEF, Glycan, Carbohydrates, Ions

Open discussion session – based on customer issues and experiences

### Day2 8.30AM-4.00PM

32 Karat software training:

User management

Training on 32 Karat depending on customer requests

Maintenance Course - Instrument Housekeeping and Troubleshooting

general housekeeping (Tips and Tricks)

electrode cleaning and removal

Interface block

detector/fibre-optic cleaning and removal

Laser old & new / Inter connect module / LIF Plug

building of Capillary cartridges

Coolant

performance check of the instrument

instrument calibration

If there are any topics or questions which you would like us to cover please forward those to <mailto:cesupport@sciex.com>

*Lunch will be provided. Please let us know if you have a dietary requirement.*