

Certificate of Analysis

Product Identification

Name: Standards Chemical Kit with Higher Concentration PPGs

Part number: 4412399

Description:

BOTTLE PART NUMBER	LOT NUMBER	COMPONENT	CONCENTRATION (+/- $\leq 2\%$)
4405227	A3322	POS PPG 1E-4M	100 μM
4405229	A3352	POS PPG 1E-5M	10 μM
4405234	A3339	NEG PPG 3E-4M	300 μM
4405236	A3337	RESERPINE 1 μM	1 μM
4405237	A3344	RESERPINE 0.167 μM	0.167 μM
4405239	A3350	RENIN 10 pmol/ μl	10 pmol/ μl
4405241	A3352	TAUROCHOLIC ACID 2 ng/ μL	2 ng/ μL

Lot Specific Information

Lot number: A3357

Manufacture date: 23 Dec 2013

Expiration date: October 2014

Quality Assurance Statement:

These standards have been manufactured using ISO9001 procedures. Concentrations are calculated gravimetrically, and balances are calibrated to NIST traceable standards in accordance with our ISO procedures. Purity from interferences is verified through mass spectrometric analysis. Solvents are minimum HPLC grade.

The information on this certificate as been reviewed and accurately reflects the manufacturing and analysis data for the product lot specified above.

Quality Assurance Associate: Brandon Munyon

Date: 23 Dec 2013

Support

Technical Support call 1 800-952-4716 or 877-740-2129 or E-mail: Support@sciex.com

This certificate may be obtained at: <http://www.sciex.com/mk/get/COFA>