

成功升级和成功大师课程所需耗材

SCIEX 培训课程遵循经验证的间隔学习法*，以提高学习效果和记忆。该培训过程在用户现场进行，是讲师指导、实验室上机实操和自主在线学习的独特结合。完成本课程后，学习者将能够熟练地调谐和校准系统、创建和优化采集方法、执行数据采集和数据处理以及对系统进行维护和故障诊断。

为了实现更好的培训效果，需要在培训前安装和配置好 SCIEX LC-MS/MS 系统。表 1 和 表 2 列出了培训期间需要使用的耗材。

表 1: 适用于所有系统的培训期间所需耗材

描述	系统	产品编码(Part number)	规格	数量
三嗪类农残混标 (100 µg/mL)	All	4376887	1 mL 安瓿瓶	1
飞诺美 Synergi 4 µm Fusion-RP 80Å HPLC 色谱柱	All	4376878	50 mm × 2.0 mm	1
注射器和注射针头	All	N/A	1 mL	1
注射器转接头	All	1008236	N/A	1
自动进样器玻璃小瓶	All	N/A	2 mL	1
玻璃瓶	All	N/A	250 mL	1
玻璃瓶	All	N/A	1 L	2
甲酸	All	N/A	N/A	N/A
LC-MS 级水	All	N/A	100 mL	N/A
LC-MS 级甲醇	All	N/A	2 L	N/A
LC-MS 级异丙醇	All	N/A	2 L	N/A
PEEK 管(红)	All	4425163	0.13 × 3000 mm	1
PEEK 管切刀	All	011281	N/A	1
尖底离心管	All	N/A	50 mL	1
Fitting PEEK 三通 0.020 in BORE	All	1006550	N/A	1

*间隔学习法是Quinnovation公司的Clark Quinn博士提出的一种学习方法。该方法理论表明，随着时间的推移，与短期内的高强度培训相比，学习并不断强化的模式能产生更好的参与度和更多的记忆保留。

表 2: 适用于特定系统的培训期间所需耗材

描述	系统	产品编码(Part number)	规格	数量
SWATH 采集性能套件	TripleTOF 系统	5045757	N/A	1
校准液套装 SCIEX X500 系统	ZenoTOF 7600 和 X500 QTOF 系统	5043301	N/A	1
LC/MS 标准品套装, 包含 High/Low PPG	所有三重四极杆和 QTRAP 系统 不包括 SCIEX 7500 系统	4406127	N/A	1
MS 单瓶调谐液	SCIEX 7500 系统	5077207	25 mL	1
MS 清洗溶剂	SCIEX 7500 系统	5077210	25 mL	1
去离子水	Echo MS 系统	N/A	1 L	N/A
Small molecule Quantitation LC-MS Performance Mixture (solution only)	Echo MS 系统	5056817	N/A	1
Echo 认证微孔板	Echo MS 系统	5077103	384 孔	1
无头 1/32 OD 红色 PEEK 管	Echo MS 系统	5075432	10	10

The SCIEX clinical diagnostic portfolio is For In Vitro Diagnostic Use. Rx Only. Product(s) not available in all countries. For information on availability, please contact your local sales representative or refer to www.sciex.com/diagnostics. All other products are For Research Use Only. Not for use in Diagnostic Procedures.

Trademarks and/or registered trademarks mentioned herein, including associated logos, are the property of AB Sciex Pte. Ltd. or their respective owners in the United States and/or certain other countries (see www.sciex.com/trademarks). Echo and Echo MS are trademarks or registered trademarks of Labcyte, Inc. in the United States and other countries, and are being used under license.

© 2023 DH Tech. Dev. Pte. Ltd. GEN-CST-02-13869-ZH-A



Headquarters
500 Old Connecticut Path | Framingham, MA 01701 USA
Phone 508-383-7700
sciex.com

International Sales
For our office locations please call the division
headquarters or refer to our website at
sciex.com/offices