

Required consumables for Success Solution for Echo[®] MS+ system with SCIEX Triple Quad 6500+ system at customer site

Our training courses follow the proven spaced learning approach to maximize learning retention. This training course includes a unique blend of self-paced eLearning, instructor led and hands-on training provided at the customer site. Upon completion of the course, you should be comfortable with creating SCIEX OS methods for the Echo[®] MS+ system with the SCIEX Triple Quad 6500+ system, handling data processing, optimizing the system and performing basic maintenance for the Echo[®] MS+ system.

For this training to be successful, the Echo[®] MS+ system must be installed and configured before the training. **Tables 1** and **2** list the consumables that must be available for use during the training. **Table 1** lists consumables that are supplied with the purchase of a new system and **Table 2** lists the consumables that you are required to provide.

Table 1: Required consumables supplied with the purchase of a new system

Description	Part number	Size	Quantity
Triazine pesticide standard mixture (100 µg/mL)	4376887	1 mL ampoule	1
Echo [®] MS Qualified microplate	5077103	384 wells	1
Fitting headless 1/32 OD tubing PEEK red	5075432	10	10

Table 2: Required consumables supplied by the customer

Description	Part number	Size	Quantity
Conical centrifuge tube	N/A	50 mL	1
Deionized water	N/A	1 L	N/A
Glass vial	N/A	2 mL	1
Glass bottle	N/A	250 mL	1
LC-MS-grade water	N/A	100 mL	N/A
LC-MS-grade methanol	N/A	2 L	N/A

The SCIEX clinical diagnostic portfolio is For In Vitro Diagnostic Use. Rx Only. Product(s) not available in all countries. For information on availability, please contact your local sales representative or refer to www.sciex.com/diagnostics. All other products are For Research Use Only. Not for use in Diagnostic Procedures.

Echo and Echo MS are trademarks or registered trademarks of Labcyte, Inc. in the United States and other countries, and are being used under license.

Trademarks and/or registered trademarks mentioned herein, including associated logos, are the property of AB Sciex Pte. Ltd. or their respective owners in the United States and/or certain other countries (see www.sciex.com/trademarks).

© 2024 DH Tech. Dev. Pte. Ltd. MKT-22719-A



Headquarters
500 Old Connecticut Path | Framingham, MA 01701 USA
Phone 508-383-7700
sciex.com

International Sales
For our office locations please call the division
headquarters or refer to our website at
sciex.com/offices