

Special trade-in promotion!



Fast and practical high-resolution separation and quantification of charged and polar analytes

Get started with the P/ACE™ MDQ Plus Capillary Electrophoresis System, with the diode array detector included, at the special price of **€49,995 + VAT**

The P/ACE™ MDQ Plus System is an automated, programmable capillary electrophoresis (CE) system designed to perform fast separations from complex samples. Separations occur in a capillary that is housed in an interchangeable cartridge with circulating liquid coolant for efficient temperature control. This finely regulated temperature control allows the use of high ionic strength buffers and large-bore capillaries, which are necessary for analyzing difficult samples.

Benefits:

- With CE, methods are faster to develop, easier to validate and less expensive to run
- Green technology reduces solvent consumption by over 100 times compared to high-performance liquid chromatography (HPLC)

Note that different detectors are available to optimize the performance of the P/ACE MDQ Plus System for your application. This special offer includes the diode array detector.

Applications include:

- Ion analysis
- Drug analysis
- Chiral isomer separation
- DNA/RNA analysis
- Protein analysis
- Oligosaccharide and sugar analysis
- Food and beverage testing

Request more information or email walter.belloni@sciex.com to get a price quote.

*Terms and conditions:

This promotion applies to the following product: B52521PDA (P/ACE™ MDQ Plus Capillary Electrophoresis System with the diode array detector). All prices include applicable customs clearance and transportation but exclude taxes. This promotion cannot be combined with other promotions or discounts and is not retroactive. No substitutions or swap-outs are allowed. The discount is not convertible to cash. All prices exclude local/state taxes. This offer cannot be used in conjunction with any other offer and has no cash value. The purchase is subject to SCIEX standard terms and conditions of sale. Other restrictions may apply. Purchase order, along with payment, must be received by 31 December 2021. Void where prohibited. Warranty: 12 months from the date of installation or 15 months from the date of shipment, whichever is earlier. Promotion available in Spain, Italy and Portugal but limitations apply. Contact SCIEX for specific availability.

The SCIEX clinical diagnostic portfolio is For In Vitro Diagnostic Use. Rx Only. Product(s) not available in all countries. For information on availability, please contact your local sales representative or refer to <https://sciex.com/diagnostics>. All other products are For Research Use Only. Not for use in Diagnostic Procedures. Trademarks and/or registered trademarks mentioned herein, including associated logos, are the property of AB Sciex Pte. Ltd. or their respective owners in the United States and/or certain other countries. © 2021 DH Tech. Dev. Pte. Ltd.

RUO-MKT-07-12500-B


The Power of Precision