

Capillary Electrophoresis Consumables

Certificate of Conformity

These products were manufactured in accordance with ISO 9001:2015.

Product Identification

Name: Low pH Phosphate SDS Sample Buffer

Kit/Reorder number: C57805

Storage: Store at 18°C to 30°C

Lot Specific Information

Lot Number:

Manufacture Date:

Expiration Date:

Quality Control Results:

The Low pH Phosphate SDS Sample Buffer complies with the regulations set forth by China Pharmacopeia's requirement of a 40 mM Phosphate Sample Buffer, pH 6.5, 1% SDS.

Quality Assurance Statement:

These reagents were verified and released by SCIEX Quality Control Laboratories and the results above were obtained.

Quality Control:
(PRINT)

Date:

Signature: _____ 

Support

Technical Support call 1 800-952-4716 or 877-740-2129 or E-mail: Support@sciex.com

This certificate may be obtained at: <http://sciex.com/tech-regulatory>