

Certificate of Conformance

These products were manufactured in accordance with compliance for ISO 9001:2015.

Product Identification

Name: Pre-Assembled Capillary Cartridge
Storage: Ambient

Part number: A55625

Lot Number:

Manufacture Date:
Expiration Date:

Quality Assurance Statement:

The components used in this assembly conform to SCIEX print and process control specifications.

Serial Number:

- 1 1105230051
- 2 1105230024
- 3 1105230027
- 4 1105230034
- 5 1105230031
- 6 1105230043
- 7 1105230011

Quality Assurance Associate:

Signature: 

Date:

Support:

Technical Support call 877-740-2129 or E-mail: sciexnow@sciex.com

This certificate may be obtained at: <https://sciex.com/support>