



Capillary Electrophoresis

Certificate of Conformance

These products were manufactured in accordance with compliance for ISO 9001:2015.

Product Identification

Name: Pre-Assembled Capillary Cartridge
Storage: Ambient

Part number: A55625

Lot Number:

Manufacture Date:

Expiration Date:

Quality Assurance Statement:

The components used in this assembly conform to SCIEX print and process control specifications.

Serial Number:

- | | | | |
|----|------------|----|------------|
| 1 | 1805230001 | 31 | 1805230061 |
| 2 | 1805230004 | 32 | 1805230063 |
| 3 | 1805230016 | 33 | 1805230065 |
| 4 | 1805230022 | 34 | 1805230071 |
| 5 | 1805230025 | 35 | 1805230073 |
| 6 | 1805230026 | 36 | 1805230075 |
| 7 | 1805230027 | 37 | 1805230081 |
| 8 | 1805230028 | 38 | 1805230082 |
| 9 | 1805230029 | 39 | 1805230085 |
| 10 | 1805230030 | 40 | 1805230086 |
| 11 | 1805230041 | 41 | 1805230088 |
| 12 | 1805230042 | 42 | 1805230090 |
| 13 | 1805230043 | 43 | 1805230097 |
| 14 | 1805230044 | 44 | 1805230098 |
| 15 | 1805230045 | 45 | 1805230102 |
| 16 | 1805230046 | 46 | 1805230104 |
| 17 | 1805230048 | | |
| 18 | 1805230049 | | |
| 19 | 1805230050 | | |
| 20 | 1805230051 | | |
| 21 | 1805230052 | | |
| 22 | 1805230053 | | |
| 23 | 1805230054 | | |
| 24 | 1805230056 | | |
| 25 | 1805230057 | | |
| 26 | 1805230058 | | |
| 27 | 1805230059 | | |
| 28 | 1805230060 | | |
| 29 | 1805230062 | | |
| 30 | 1805230064 | | |

Quality Assurance Associate:

Signature: *sk*

Date:

Support:

Technical Support call 877-740-2129 or E-mail: sciexnow@sciex.com

This certificate may be obtained at: <https://sciex.com/support>