



Certificate of Conformance

These products were manufactured in accordance with compliance for ISO 9001:2015.

Product Identification

Name: Pre-Assembled Capillary Cartridge
Storage: Ambient

Part number: A55625

Lot Number:


Manufacture Date:
Expiration Date:

Quality Assurance Statement:

The components used in this assembly conform to SCIEX print and process control specifications.

Serial Number:

- 1 2311230029
- 2 2311230030
- 3 2311230047

Quality Assurance Associate:
Signature: 

Date:

Support:

Technical Support call 877-740-2129 or E-mail: sciexnow@sciex.com

This certificate may be obtained at: <https://sciex.com/support>